

System specifications

Symbia Pro.specta **SPECT/CT**

siemens-healthineers.com/SymbiaPro.specta

Symbia Pro.specta SPECT/CT

System specifications

System hardware (standard)
Open gantry design (70 cm/27.6 in)
Tin Filter
Stellar CT detector
Two high-definition digital SPECT detectors
Low-profile 3/8" or 5/8" detectors
Detector configuration, including 180°, 90°, 76°, gurney, standing, sitting, out-facing
Automatic body contouring
Detector touchpad sensors
Lightrail sensors and lightrail touchpads
Circular and non-circular body-contour orbit
Touch screen gantry display (Scan&GO) used for the following functions: <ul style="list-style-type: none"> – Patient name and date of birth – Patient positioning with window and persistence adjustment – Acquisition parameter display (time remaining, count rate, etc.) – Camera information (detector and bed positions) – Gantry control (reconfiguration, collimator change, offset zoom, reset, etc.) – Unload patient
Patient bed with 227 kg (500 lb) capability

System hardware (standard)
Patient bed pivot for rail-free access of sitting/standing patients, wheelchairs, imaging tables, gurneys, and hospital beds
Integrated calibration source holder
Rear bed with pallet flex prevention
Carbon fiber pallet
Patient comfort accessories (SPECT/CT head holder and cushion set, head arm support, head rest, patient restraint set, body wraps set, whole-body armrest set)
Acquisition workplace with multilingual graphical user interface, monitor, keyboard, and mouse plus full DICOM archiving, and external USB 3.0 disks support
Intuitive hand controller with easy-to-use descriptive controls
24-in/58-cm flat screen monitor
Operator manuals
CT quality assurance phantom

Symbia Pro.specta SPECT/CT

System specifications

System hardware (optional)
Caudal and cephalic tilt
Low-energy high-resolution (LEHR) collimator
Low-energy all-purpose (LEAP) collimator
Low-energy ultra-high-resolution (LEUHR) collimator
Medium-energy low-penetration (MELP) collimator
Low-penetration high-resolution (LPHR) collimator
High-energy (HE) collimator
Pinhole collimator (4, 6, or 8 mm aperture)
IQ•SPECT™ with SMARTZOOM™ collimator
SMARTZOOM HRX Collimator
Integrated Collimator Changer (ICC)
Automatic Collimator Changer (ACC)
Collimator cart(s)
Automatic Quality Control (AQC)
AutoQC source kit (base or premium)
NIST-traceable precision calibration source (⁵⁷ Co, ⁷⁵ Se, ¹¹³ Sn)
Internal ECG (AHA or IEC)
External ECG (AHA or IEC)
Patient handling system (PHS) extended pivot
Dual monitor
Plan&GO

System hardware (optional)
Radiation therapy pallet
Mammography pallet
Pediatric pallet
Under- or over-floor PHS cable
Seismic installation kit (single or dual anchor)
Integrated Electronics Cabinet (IEC)
Uninterruptible Power Supply (UPS)
Power bundle (IEC and UPS)
High-performance integrated computer system (ICS)
FAST IRS (image reconstruction system) ¹
Ultra-FAST IRS ²
4-quadrant phantom
Rectangular fillable flood source
Deluxe SPECT phantom
American College of Radiology (ACR) SPECT phantom kit
Computer workstation table
Swivel tilt chair
Additional user manuals
10-mCi rectangular source
CT high-speed 0.33s ³

System specifications

System software (standard)
SPECTsyngo software
Integrated acquisition and reconstruction
Planar (static) acquisition
Dynamic acquisition
Whole-body acquisition
Whole-body SPECT acquisition
Single-bed SPECT acquisition
Gated acquisition
Gated SPECT acquisition
Dynamic SPECT acquisition
Filtered back projection
Flash 3D iterative reconstruction imaging
Scatter correction
Reconstruction with CT frame-of-reference
CT attenuation correction
Smart Remote Services (SRS)
Gantry display and interaction at acquisition workplace
Cardiac retrospective gating
Automatic intra-recon motion correction
Topogram, spiral, sequential acquisition modes
Extended field of view (FOV) (HD FoV)
SAFIRE
IVR (Interleaved Volume Reconstruction)
CARE kV ³
CARE Child
CARE Filter
CARE Bolus CT
CARE Topo
CARE Dose4D™
X-CARE

System software (standard)
FAST Planning
FAST ROI
FAST kV
SureView™
Flexible kV (80, 110, 130 kV) (130 kV not with CT Intro functionality)
Flexible kV (70-140 kV in ranges of 10 kV) ³
CT DICOM structured dose report
syngo archiving and network
Dynamic planar real-time activity curve
Video capture and editing tool
Exam designer
WorkStream4D™ (direct 3D-reconstruction)
Filming
myExam Companion (includes myExam Compass)
Check&GO
Recon&GO
MI View&GO
syngo® Security Package
Endoscopic view
Flex Dose Profile

System specifications

System software (optional)
xSPECT reconstruction with CT frame-of-reference ²
xSPECT Bone™ with zone map (tissue classification) ²
xSPECT Quant™ ^{99m} Tc ²
xSPECT Quant ¹⁷⁷ Lu ²
xSPECT Quant ¹²³ I ²
xSPECT Quant ¹¹¹ In ²
xSPECT Quant ¹³¹ I ²
Broad Quantification™
Dose calibrator cross-calibration capability for unbiased SUV quantification
TrueCalc™ high-count-rate detector technology
3D measured collimator characterization for IQ•SPECT (hole, shape, and size)
AutoPlanar
NM Basic cardiac package (respiratory motion correction, continuous gated acq, gated AC/SC)
SMART Neuro AC
Organ Processing for Symbia™
Planar half-time imaging
syngo Security Package
syngo.MI Cardiology 4DM
syngo.CT Extension Corridor4DM™
syngo.SPECT Cardiology Cedars
syngo.CT Extension Cedars
syngo.MI Cedars Reporting
syngo.PET Cardiology Cedars

System software (optional)
iMAR
CT Lung CAD package
CT syngo Volume Perfusion Neuro
CT Neuro DSA
CT syngo Dental
CT Dual Energy Base Package
CT CaScoring package
syngo Expert-i
myExam Cockpit
CT Cardio Base Package ³
CT Cardio Advanced Package ³
CT Dual Energy Base Package
IQ•SPECT software

SPECT specifications

Gantry dimensions	
Height	223.5 cm (7 ft 4 in)
Width	235.5 cm (7 ft 8.7 in)
Depth	206.9 cm (6 ft 9.5 in)
Axis of rotation (from floor)	104.5 cm (3 ft 4.7 in)
Minimum/maximum patient opening (HE collimator)	12 cm (4.7 in)/65.4 cm (25.7 in)
Minimum/maximum patient opening (LEHR collimator)	19.2 cm (7.6 in)/72.6 cm (28.6 in)
Scan&GO	39.6 cm (15.6 in) flat panel color LCD display
Tunnel opening	70 cm aperture (27.6 in)
Tunnel length	106.37 cm (41.88 in)
Distance between SPECT and CT FOV ⁴	116.7 cm (45.9 in)

SPECT acquisition	
Energy range (photopeak window center)	35-588 keV
Acquisition modes	Static, dynamic, gated, SPECT, gated SPECT, dynamic SPECT, composite SPECT, whole-body, whole-body SPECT, SPECT/CT

Nuclear medicine and SPECT acquisition parameters	Static
Time	50-32,000,000 ms
Counts	1-2,000,000 cts
Zoom	1.00, 1.231, 1.455, 1.778, 2.00, 2.286, 2.667, 3.20
Matrix	64 x 64, 128 x 128, 256 x 256, 512 x 512, 1,024 x 1,024
Patient position	Head first supine (HFS), head first prone (HFP), feet first supine (FFS), feet first prone (FFP), upright, head left supine (HLS), head right supine (HRS), head left prone (HLP), head right prone (HRP)
Detectors	Detector 1, detector 2, both
Detector configuration	180°, 90°, 76°, gurney, outer room left, outer room right
Allowable collimators ⁵	LEHR, LPHR, LEAP, LEUHR, ME, HE, pinhole

SPECT specifications

Nuclear medicine and SPECT acquisition parameters	
	Dynamic
Time	1-1,440 min
Time per frame per phase	500-32,000,000 ms (dynamic planar); 10-32,000,000 msec (cardiac dynamic planar)
Zoom	1.00, 1.231, 1.455, 1.778, 2.00, 2.286, 2.667, 3.20
Matrix	64 x 64, 128 x 128, 256 x 256
Number of phases	1-8 phases
Patient position	HFS, HFP, FFS, FFP, upright, HLS, HRS, HLP, HRP
Detectors	Detector 1, detector 2, both
Detector configuration	180°, 90°, 76°, gurney, outer room left, outer room right
Acquire with R-wave gate	Selectable
Acquire with statics	Selectable
Allowable collimators ⁵	LEHR, LPHR, LEAP, LEUHR, MELP, HE, pinhole
	Gated
Time	500-32,000,000 ms
Counts	1-2,000,000 cts
Zoom	1.00, 1.231, 1.455, 1.778, 2.00, 2.286, 2.667, 3.20
Matrix	64 x 64, 128 x 128
Number of gates	6-32
Patient position	HFS, HFP, FFS, FFP, upright, HLS, HRS, HLP, HRP
Detectors	Detector 2
Detector configuration	180°, 90°, 76°, gurney, outer room left, outer room right
Heart beat framing	100% forward (PCNT)
Reject PVC beats	Selectable
Beats to reject post PVC	0-6
PVC threshold (bpm)	0-40% of beat upper limit
Allowable collimators ⁵	LEHR, LPHR, LEAP, LEUHR, MELP, HE, pinhole

SPECT specifications

Nuclear medicine and SPECT acquisition parameters	Whole-body
Scan speed with autocontour	30-600 mm/min
Zoom	1.00
Matrix	256 x 512, 256 x 1,024, 512 x 1,024
Scan length	387-2,000 mm
Patient position	FFS, FFP
Detectors	Detector 1, detector 2, both
Detector configuration	180°
Autocontour	Selectable
Allowable collimators ⁵	LEHR, LPHR, LEAP, LEUHR, MELP, HE
	SPECT
Time	500-32,000,000 ms
First view by counts	1-2,000,000 kcts
Zoom	613.8 mm/CT width (with AC CT); 1.00 (without AC CT)
Matrix	64 x 64, 128 x 128, 256 x 256
Maximum number of views	360 per head
Patient position	HFS, HFP, FFS, FFP
Detectors	Detector 1, detector 2, both
Detector configuration	180°, 90°, 76°, IQ•SPECT
Orbit	Circular (180°, 90°), NCO (180°, 90°, 76°), NCO-prescan (90°), cardio-centric
Start angle	-179°-180°
Mode	Step and shoot, continuous
Degrees of rotation	90° (only 90°), 104° (76°, IQ•SPECT), 180° (only 180°), 360° (90° and 180°)
Rotation direction	Clockwise, counterclockwise
Allowable collimators ⁵	LEHR, LPHR, LEAP, LEUHR, MELP, HE

SPECT specifications

Nuclear medicine and SPECT acquisition parameters	Composite SPECT
Number of rotations	1-45
Rotation time	10-3,600 sec (continuous); 8-11,5200,000 (step & shoot)
Zoom	613.8 mm/CT width (with AC CT); 1.00 (without AC CT)
Matrix	64 x 64, 128 x 128, 256 x 256
Start angle	-179°-180°
Patient position	HFS, HFP, FFS, FFP
Detectors	Both
Detector configuration	180°
Orbit	Circular, NCO
Mode	Step and shoot, continuous
Rotation direction	Clockwise, counterclockwise
Degrees of rotation	180°, 360°
Allowable collimators ⁵	LEAP, LPHR, LEHR, LEUHR, MELP, HE
	Dynamic SPECT
Number of rotations	1-45
Rotation time	10-3,600 sec (continuous); 8-11,5200,000 (step & shoot)
Number of phases	1–total number or rotations specified
Recon images	1–total number or rotations per phase
Zoom	613.8 mm/CT width (with AC CT); 1.00 (without AC CT)
Matrix	64 x 64, 128 x 128, 256 x 256
Start angle	-179°-180°
Patient position	HFS, HFP, FFS, FFP
Detectors	Both
Detector configuration	180°, 90°, 76°, IQ•SPECT
Orbit (at detector configuration)	Circular (180°, 90°), NCO (180°, 90°, 76°), NCO-prescan (90°), cardio-centric
Mode	Step and shoot, continuous
Rotation direction	Clockwise, counterclockwise
Degrees of rotation (at detector configuration)	90° (only 90°), 140° (76°, IQ•SPECT), 180° (only 180°), 360° (90° and 180°)
Allowable collimators ⁵	LEAP, LPHR, LEHR, LEUHR, MELP, HE

SPECT specifications

Nuclear medicine and SPECT acquisition parameters	Gated SPECT
Time	500-32,000,000 ms
Zoom	613.8 mm/CT width (with AC CT); 1.00 (without AC CT)
Matrix	64 x 64, 128 x 128, 256 x 256
Number of gates	6-32
Patient position	FFS, FFP
Detectors	Detector 1, detector 2, both
Detector configuration	180°, 90°, 76°, IQ•SPECT
Orbit	Circular (180°, 90°), NCO (180°, 90°, 76°), NCO-prescan (90°), cardio-centric
Start angle	-179°-180°
Mode	Step and shoot, continuous
Degrees of rotation	90° (only 90°), 104° (76°, IQ•SPECT), 180° (only 180°), 360° (90° and 180°)
Rotation direction	Clockwise, counterclockwise
Heart beat framing	100% forward (PCNT)
Reject PVC beats	Selectable
Beats to reject post PVC	0-6
PVC threshold (bpm)	0-40% of beat upper limit
Allowable collimators ⁵	LEHR, LPHR, LEAP, LEUHR, MELP, HE

SPECT specifications

SPECT motions	
Average autocontour distance	1.1 cm (0.45 in)
Maximum radial and lateral speed	72 cm/min (28.3 in/min)
Maximum lateral position left/right	37.5 cm (14.7 in)/10 cm (4 in)
Maximum clockwise/counterclockwise rotation detector 1	405°/-135°
Ring rotation range	540° (from -135° to +405°)
Rotational uniformity	Yes
Rotational accuracy	0.1°
Rotational speed	0.03-3.0 RPM
Center of rotation	≤0.25 pixel (64 x 64 matrix)
Maximum caudal tilt	+16°/-16°
Patient bed	
Width	68.6 cm (27.0 in)
Length	243.8 cm (8 ft 1 in)
Weight without ICC/ACC	860.4 kg (1,896.8 lb)
Height	101.5 cm (3 ft 4 in)
Vertical motion range	48.0-109.0 cm (19-43 in)
Vertical speed	4-12 mm/sec (.16-.47 in/sec)
Pallet material	Carbon fiber
Pallet thickness	15 mm (.6 in)
Pallet width	39.8 cm (15.7 in)
Attenuation at 140 keV	<10%, <12% (with accessories)
Maximum patient weight	227 kg (500 lb)
Maximum deflection of patient pallet	<2.0 mm (<0.08 in) for 90.7 kg (200 lb ± 100) patient
Maximum scan length in whole-body mode	200 cm (6 ft 6.7 in)
Horizontal motion accuracy	+/- 0.5 mm (0.02 in)
Minimum/maximum horizontal speed	0.5-100 mm/sec (0.02-3.94 in/sec)

SPECT specifications

Optional pallets		
Pediatric	Material	Carbon fiber composite
	Thickness	0.5 cm (0.18 in)
	Width	25.4 cm (<10 in)
	Length	74.9 cm (29.5 in)
	Weight	7.0 kg (15.4 lb)
	Attenuation at 140 keV	<10%
	Maximum patient weight	27 kg (60 lb)
Scintimammography	Material	Carbon fiber composite
	Thickness	1.6 cm (0.63 in)
	Width	35.6 cm (14 in)
	Length	183.5 cm (72.3 in)
	Weight	7.6 kg (16.8 lb)
	Attenuation at 140 keV	<10%
	Maximum patient weight	135 kg (300 lb)
Radiotherapy planning	Material	Carbon fiber composite
	Thickness	1.5 cm (0.6 in)
	Width	53 cm (20.9 in)
	Length	203.5 cm (80.1 in)
	Weight	9.0 kg (19.8 lb)
	Attenuation at 140 keV	<10%
	Maximum patient weight	227 kg (500 lb)
Rear pallet support		
Width	26.9 cm (10.6 in)	
Length	113.7 cm (3 ft 8.8 in)	
Weight	182.9 kg (403.3 lb)	

SPECT specifications

ECG trigger	
Integration	Internal (inside patient bed) or external
Framing modes	100% forward (PCNT)
Bad beat rejection	Yes
Criteria for framing images	Frames/R-R interval
Beat acceptance window	Automatic or manual selection

Collimator cart	
Height	1,016 mm (40.0 in)
Width	822 mm (32.3 in)
Depth	1,091 mm (43.0 in)
Weight (without collimators)	173.8 kg (383.2 lb)

Detector dimensions	
FOV	53.3 x 38.7 cm (21 x 15.25 in)
Diagonal FOV	65.9 cm (25.9 in)

Crystal	
Size	59.1 x 44.5 cm (23.25 x 17.5 in)
Diagonal	73.9 cm (29.1 in)
Thickness	9.5 mm (3/8 in) or 15.9 mm (5/8 in)

Photomultiplier tubes	
Total number	59
Diameter	53-7.6 cm (3 in) and 6-5.1 cm (2 in)
Type	Bialkali high-efficiency box-type dynodes
Array	Hexagonal
Sampling rate	30.0 MHz

Detector shielding	
Back	Typical 15.88 mm (0.625 in) Minimum 9.5 mm (0.374 in)
Sides	12.7 mm (0.5 in)
Minimum/maximum in patient direction ⁶	27.9/36.4 mm (1.1/1.435 in)
Brain reach ⁷	7.6 cm (3 in)

SPECT specifications

Detector ⁸	3/8" (0.9525 cm)	5/8" (1.5875 cm)
Intrinsic spatial resolution		
Full width at half maximum (FWHM) in central field of view (CFOV)	≤3.84 mm	≤4.54 mm
FWHM in useful field of view (UFOV)	≤3.94 mm	≤4.64 mm
Full width at tenth maximum (FWTM) in CFOV	≤7.54 mm	≤8.70 mm
FWTM in UFOV	≤7.74 mm	≤8.90 mm
Intrinsic spatial linearity		
Differential in CFOV	≤0.24 mm	≤0.24 mm
Differential in UFOV	≤0.24 mm	≤0.24 mm
Absolute in CFOV	≤0.44 mm	≤0.54 mm
Absolute in UFOV	≤0.69 mm	≤0.99 mm
Intrinsic energy resolution		
FWHM in UFOV	≤9.85%	≤9.85%
Intrinsic flood field uniformity (uncorrected)		
Differential in CFOV	≤2.54%	≤2.54%
Differential in UFOV	≤2.74%	≤2.74%
Integral in CFOV	≤2.94%	≤2.94%
Integral in UFOV	≤3.74%	≤3.74%
Intrinsic count rate performance in air		
Maximum count rate	≥460 kcps	
Maximum count rate (@15% window)	≥310 kcps	≥310 kcps
Intrinsic spatial resolution at 75 kcps		
FWHM in UFOV	≤4.1 mm	≤4.64 mm
FWTM in UFOV	≤7.84 mm	≤8.9 mm
Intrinsic flood field uniformity at 75 kcps (uncorrected)		
Differential in CFOV	≤2.54%	≤2.54%
Differential in UFOV	≤2.74%	≤2.74%
Integral in CFOV	≤2.94%	≤2.94%
Integral in UFOV	≤3.74%	≤3.74%
Multiple window spatial registration	≤0.64 mm	≤0.99 mm

SPECT specifications

High-count-rate performance ⁹	3/8" (0.9525 cm)	5/8" (1.5875 cm)
Detector specifications at 310 kcps¹⁰		
Intrinsic flood field uniformity (uncorrected)		
Differential in CFOV	≤3.04%	–
Differential in UFOV	≤3.24%	–
Integral in CFOV	≤3.44%	–
Integral in UFOV	≤4.24%	–
Intrinsic energy resolution		
FWHM in CFOV	≤11.85%	–
Stability of energy peak position		
Change of peak position (70 kcps versus 30 kcps) ¹⁰	≤± 0.5%	–
System spatial resolution without scatter (LEHR at 10 cm)		
FWHM in CFOV	≤8.04 mm	–
FWTM in CFOV	≤14.64 mm	–
Detector with collimator⁸		
System spatial resolution without scatter (LEHR at 10 cm)		
FWHM in CFOV	≤7.54 mm	≤7.84 mm
FWTM in CFOV	≤13.64 mm	≤14.94 mm
System spatial resolution with scatter (LEHR at 10 cm)		
FWHM in CFOV	≤8.34 mm	≤8.9 mm
FWTM in CFOV	≤18.64 mm	≤19.5 mm
System planar sensitivity (LEHR at 10 cm)		
Absolute ^{99m} Tc	202 cpm/μCi +/- 10%	225 cpm/μCi +/- 10%
System planar sensitivity (MELP at 10 cm)		
Absolute ¹¹¹ In	430 cpm/μCi +/- 20%	565 cpm/μCi +/- 20%

SPECT specifications

Detector with collimator tomographic ⁸	3/8" (0.9525 cm)	5/8" (1.5875 cm)
Reconstructed spatial resolution without scatter correction at 15-cm radius (LEHR)		
	Filtered back projection	
Central transaxial	≤10.24 mm	–
Central axial	≤10.84 mm	–
Peripheral radial	≤9.84 mm	–
Peripheral tangential	≤8.44 mm	–
Peripheral axial	≤9.04 mm	–
Reconstructed spatial resolution without scatter correction at 15-cm radius (LEHR)		
	Flash 3D iterative reconstruction	
Central transaxial	≤4.44 mm	–
Central axial	≤4.44 mm	–
Peripheral radial	≤4.04 mm	–
Peripheral tangential	≤3.94 mm	–
Peripheral axial	≤4.24 mm	–
Reconstructed spatial resolution with scatter (LEHR)		
	Filtered back projection	
Center	≤10.74 mm	≤11.54 mm
Radial	≤10.94 mm	≤12.04 mm
Tangential	≤7.94 mm	≤8.84 mm
Reconstructed spatial resolution with scatter (LEHR)		
	Flash 3D iterative reconstruction	
Center	≤5.84 mm	–
Radial	≤5.04 mm	–
Tangential	≤4.14 mm	–
Average volume sensitivity per axial centimeter		
LEHR, ^{99m} Tc	12,000 (cts/sec)/(MBq/cm ²)	–
Detector-to-detector sensitivity variation		
LEHR, ^{99m} Tc	≤5.0%	–
Detector with collimator whole-body scanning⁸		
	3/8" (0.9525 cm)	5/8" (1.5875 cm)
Whole-body system spatial resolution without scatter at 10-cm/min scan speed (LEHR at 10 cm)		
FWHM perpendicular	≤7.54 mm	–
FWHM parallel	≤7.94 mm	–
FWTM perpendicular	≤14.04 mm	–
FWTM parallel	≤14.24 mm	–

SPECT specifications

Collimators	LEHR	LPHR	LEAP	LEUHR	MELP	HE	SMARTZOOM	SMARTZOOM HRX
	Low Energy High Resolution	Low Penetration High Resolution	Low Energy All Purpose	Low Energy Ultra-high Resolution	Medium Energy Low Penetration	High Energy	IQ•SPECT	Ultra-high Resolution
Isotope	^{99m} Tc	¹²³ I	^{99m} Tc	^{99m} Tc	⁶⁷ Ga	¹³¹ I	^{99m} Tc	^{99m} Tc, ¹²³ I
Hole shape	Hex	Hex	Hex	Hex	Hex	Hex	Hex	Hex
Number of holes (x1,000)	148	86	90	146	14	8	48	65
Hole length	24.05 mm	35.0 mm	24.05 mm	35.8 mm	40.64 mm	59.7 mm	40.25 mm	40.25 mm
Septal thickness	0.16 mm	0.2 mm	0.2 mm	0.13 mm	1.14 mm	2 mm	0.2-0.4 mm	0.26 mm
Hole diameter across the flats	1.11 mm	1.5 mm	1.45 mm	1.16 mm	2.94 mm	4 mm	1.9 mm	1.5 mm
Sensitivity at 10 cm ¹¹	202 cpm/μCi	251 cpm/μCi (167 cpm/μCi (^{99m} Tc))	330 cpm/μCi	88 cpm/μCi	275 cpm/μCi	135 cpm/μCi	285 cpm/μCi	151 cpm/μCi (^{99m} Tc), 180 cpm/μCi (¹²³ I)
Sensitivity at 24 cm ¹²								295 cpm/μCi (^{99m} Tc), 301 cpm/μCi (¹²³ I) at 24 cm
Sensitivity at 28 cm ¹²								371 cpm/μCi (^{99m} Tc), 369 cpm/μCi (¹²³ I) at 28 cm
Geometric resolution at 10 cm	6.4 mm	6.4 mm	8.3 mm	4.6 mm	10.8 mm	13.2 mm	6.95 mm	5.41 mm
System resolution at 10 cm	7.5 mm	8.0 mm	9.4 mm	6.0 mm	12.5 mm	13.4 mm	7.4 mm ¹³	6.1 mm (^{99m} Tc), 6.4 mm (¹²³ I)
Calculated penetration	1.5%	1.2%	1.9%	0.8%	1.2%	3.5%	N/A	N/A

SPECT specifications

Pinhole collimator ¹⁴	Isotope		
	^{99m} Tc	¹²³ I	¹³¹ I
Hole shape	Round	Round	Round
Number of holes	1	1	1
Cone aperture	4 mm	4 mm	4 mm
	6 mm	6 mm	6 mm
	8 mm	8 mm	8 mm
Cone length	209.96 mm	209.96 mm	209.96 mm
Diameter at base of cone (approximate)	220.14 mm	220.14 mm	220.14 mm
Sensitivity at 10 cm with 4 mm	123 cpm/μCi	111 cpm/μCi	67 cpm/μCi
Sensitivity at 10 cm with 6 mm	271 cpm/μCi	243 cpm/μCi	133 cpm/μCi
Sensitivity at 10 cm with 8 mm	478 cpm/μCi	426 cpm/μCi	221 cpm/μCi
Geometric resolution at 10 cm with 4 mm	6.2 mm	6.3 mm	7.5 mm
Geometric resolution at 10 cm with 6 mm	9.3 mm	9.3 mm	10.6 mm
Geometric resolution at 10 cm with 8 mm	12.3 mm	12.4 mm	13.6 mm
System resolution at 10 cm with 4 mm	6.64 mm	6.64 mm	7.64 mm
System resolution at 10 cm with 6 mm	9.54 mm	9.54 mm	10.74 mm
System resolution at 10 cm with 8 mm	12.54 mm	12.54 mm	13.74 mm

CT specifications

Gantry dimensions			
	Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
Scan field	50/70 ¹⁵	50/70 ¹⁵	50/70 ¹⁵
Rotation time	0.33 s ¹⁶		
	0.5 s	0.8 s	0.8 s
	1.0 s	1.0 s	1.0 s
	1.5 s	1.5 s	1.5 s
Temporal resolution	165 ms, 83 ms (bisegment) ¹⁶	400 ms	400 ms ¹⁶
Data acquisition system			
Maximum number of slices/rotation	32 (acquired slices)	32 (acquired slices)	16 (acquired slices)
	64 (reconstructed slices)	64 (reconstructed slices)	32 (reconstructed slices)
Number of physical detector rows	32	32	16
Number of detector elements	24,576	24,576	12,288
Number of projections	1,536	1,536	1,536
Sequence acquisition modes	32 x 0.7 mm		
	Sn 32 x 0.7 mm		
	2 x 5 mm	32 x 0.7 mm	16 x 0.7 mm
	1 x 10 mm	Sn 32 x 0.7 mm	Sn 16 x 0.7 mm
	1 x 5 mm	1 x 10 mm	1 x 10 mm
	3 x 5 mm	1 x 5 mm	1 x 5 mm
Spiral acquisition modes	32 x 0.7 mm	32 x 0.7 mm	16 x 0.7 mm
	Sn 32 x 0.7 mm	Sn 32 x 0.7 mm	Sn 16 x 0.7 mm

CT specifications

Tube assembly	Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
Tube	Athlon™ X-ray tube	Chronon™ X-ray tube	Chronon X-ray tube
Tube current	13-625 mA 13-825 mA ¹⁶	13-240 mA 13-400 mA ¹⁶	13-240 mA 13-400 mA ¹⁶
Tube voltage	70-140 kV in 10 kV steps	80 kV, 110 kV, 130 kV, Sn110 kV, Sn130 kV	80 kV, 110 kV, 130 kV ¹⁶ , Sn110 kV, Sn130 kV ¹⁶
Tube anode heat storage capacity	7.0 MHU; equivalent to 17.5 MHU with SAFIRE	3.5 MHU; equivalent to 8.75 MHU with SAFIRE	3.5 MHU; equivalent to 8.75 MHU with SAFIRE
Focal spot size according to IEC 60336	0.8 x 0.8/7° 1.0 x 1.2/7°	0.8 x 0.4/8° 0.8 x 0.7/8°	0.8 x 0.4/8° 0.8 x 0.7/8°
CARE Filter			
CARE Filter tube	Equivalent to 5.5 mm Al @ 140 kV	Equivalent to 5.5 mm Al @ 140 kV	Equivalent to 5.5 mm Al @ 140 kV
CARE Filter beam limiting device	Equivalent to 0.5 mm in the isocenter	Equivalent to 0.5 mm in the isocenter	Equivalent to 0.5 mm in the isocenter
Generator			
Maximum power	75 kW; equivalent 187 kW max. generator power with SAFIRE	32 kW; equivalent 80 kW max. generator power with SAFIRE	32 kW; equivalent 80 kW max. generator power with SAFIRE
Topogram			
Length (maximum)	1,950 mm	1,950 mm	1,950 mm
Scan times	1.9 s-20.2 s	1.9 s-20.2 s	1.9 s-20.2 s
Views	a.p., p.a., lateral	a.p., p.a., lateral	a.p., p.a., lateral
Sequence acquisition			
Reconstructed slice widths	1.5, 2, 3, 4, 5, 6, 7, 8, 10 mm	1.5, 2, 3, 4, 5, 6, 7, 8, 10 mm	1.5, 2, 3, 4, 5, 6, 7, 8, 10 mm
Scan times full scan (360°)	0.33 ¹⁶ , 0.5, 1.0 s	0.8, 1.0, 1.5 s	0.8, 1.0, 1.5 s
Partial scan times (240°)	0.22 ¹⁶ , 0.33, 0.67 s	0.54, 0.67, 1.01 s	0.54, 0.67, 1.01 s

CT specifications

Multislice spiral acquisition	Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
Reconstructed slice widths	0.6, 0.8, 1, 1.5, 2, 3, 4, 5, 6, 7, 8, 10 mm	0.6, 0.8, 1, 1.5, 2, 3, 4, 5, 6, 7, 8, 10 mm	0.6, 0.8, 1, 1.5, 2, 3, 4, 5, 6, 7, 8, 10 mm
Scan times full scan (360°)	0.33 ¹⁶ , 0.5, 1.0 s	0.8, 1.0, 1.5 s	0.8, 1.0, 1.5 s
Reconstruction increment	Min. 0.1 mm	Min. 0.1 mm	Min. 0.1 mm
Pitch factor	0.03-1.5	0.09-1.5	0.09-1.5
Spiral scan time maximum	Max. 300 s	Max. 300 s	Max. 300 s
CT scan range	Max. 1,886.2 mm	Max. 1,913.7 mm	Max. 1,930.6 mm
Continuous scan length and SPECT/CT co-scan range	Max. 1,843 mm	Max. 1,843 mm	Max. 1,843 mm
Continuous scan length and SPECT/CT co-scan range minimum small room	1,483.4 mm	1,496.8 mm	1,438.6 mm
Image reconstruction			
Real-time display	512 x 512	512 x 512	512 x 512
Slice thickness	0.6-10 mm	0.6-10 mm	0.6-10 mm
Scan field	50 cm/19.7 in 70 cm/27.5 in w/HD FoV Pro ¹⁵	50 cm/19.7 in 70 cm/27.6 in w/HD FoV Pro ¹⁵	50 cm/19.7 in 70 cm/27.6 in w/HD FoV Pro ¹⁵
Reconstruction field	5-50 cm/1.9 in-19.7 in 5-70 cm/1.9 in-27.5 in w/HD FoV Pro ¹⁵	5-50 cm/1.9 in-19.69 in 5-70 cm/1.9 in-27.7 in w/HD FoV Pro ¹⁵	5-50 cm/1.9 in-19.69 in 5-70 cm/1.9 in-27.7 in w/HD FoV Pro ¹⁵
Reconstruction time	13 fps for FBP, 9 fps for IR	13 fps for FBP, 9 fps for IR	13 fps for FBP, 9 fps for IR
Reconstruction matrix	512 x 512	512 x 512	512 x 512
Hounsfield unit (HU) scale	-8,192 to +57,343	-8,192 to +57,343	-8,192 to +57,343

CT specifications

Phantom CATPHAN (16 cm)	Symbia Pro.specta X3	Symbia Pro.specta Q3
Object size	3 mm	3 mm
Contrast difference	3 HU	3 HU
Dose at surface	12.84 mGy	13.74 mGy
Technique	1.0 s, 10 mm, 130 kV	1.0 s, 10 mm, 130 kV

Phantom CATPHAN (20 cm)	Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
Object size	5 mm	5 mm	5 mm
Contrast difference	3 HU	3 HU	3 HU
Dose at surface	12.22 mGy	12.84 mGy	13.74 mGy
Technique	1.0 s, 10 mm, 120 kV	1.0 s, 10 mm, 130 kV	1.0 s, 10 mm, 130 kV

Dose, CTDI100 Values ¹⁷ (in mGy/100 mAs)	Symbia Pro.specta X7				Symbia Pro.specta X3		Symbia Pro.specta Q3	
Kv	70	80	110	140	110	130	110	130
16 cm A	4.0	6.3	15.7	27.1	13.1	19.7	13.4	20.3
16 cm B	4.5	6.8	16.4	27.9	14.2	20.9	14.9	22.3
32 cm A	0.9	1.6	4.6	8.7	3.9	6.1	4.0	6.4
32 cm B	2.2	3.5	8.9	15.9	7.8	11.9	8.3	12.8

A: at center

B: 1 cm below the surface

Technique:

- PMMA phantom
- Absorbed dose for reference material air
- Max Deviation:
 - ±20% for tube currents ≥25 mA
 - ±30% for tube currents <25 mA
- Expected deviation:
 - ±10% for tube currents ≥25 mA
 - ±20% for tube currents <25 mA

High-contrast resolution	Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
2% MTF (±10%)	15.1 lp/cm	15.1 lp/cm	15.0 lp/cm
10% MTF (±10%)	14.6 lp/cm	14.6 lp/cm	14.5 lp/cm
50% MTF (±10%)	12.0 lp/cm	12.0 lp/cm	11.8 lp/cm
Technique: Tungsten wire in air	160 mA, 120 kV, 1 s, 5 mm	105 mA, 130 kV, 1.5 s, 5 mm	105 mA, 130 kV, 1.5 s, 5 mm

xSPECT advanced specifications

Advanced bone imaging	Symbia Pro.specta X7	Symbia Pro.specta X3	
Context-based information	Yes, applied to ^{99m} Tc diphosphonate bone SPECT	Yes, applied to ^{99m} Tc diphosphonate bone SPECT	
Extra modality information	Zone map (a map with up to 6 tissue zones)	Zone map (a map with up to 6 tissue zones)	
CT zone classification	Cortical bone, spongius bone, soft tissue, air, adipose (fat), metal	Cortical bone, spongius bone, soft tissue, air, adipose (fat), metal	
Reconstruction software	xSPECT Bone	xSPECT Bone	
Reconstruction matrix size	128 x 128, 256 x 256	128 x 128, 256 x 256	
Attenuation map	Linear attenuation coefficients @ 140 keV	Linear attenuation coefficients @ 140 keV	


SPECT quantification	Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
System calibration source	NIST traceable precision ⁵⁷ Co source, ⁷⁵ Se source, ¹¹³ Sn source, point sources for specific isotope	NIST traceable precision ⁵⁷ Co source, ⁷⁵ Se source, ¹¹³ Sn source, point sources for specific isotope	Point sources for specific isotope
System calibration procedure	Every 6 months	Every 6 months	Every 6 months
Data format	Data is saved in PET format	Data is saved in PET format	Data is saved in PET format
Reconstruction software	xSPECT Quant	xSPECT Quant	
Volumetric analysis software	MI View&GO, syngo.via	MI View&GO, syngo.via	MI View&GO, syngo.via
Quantitative volumetric analysis	In units of Bq/ml, SUV or count-rate-per-voxel	In units of Bq/ml, SUV or count-rate-per-voxel	
Absolute quantification	xSPECT Quant ^{99m} Tc, ¹²³ I, ¹¹¹ In, ¹⁷⁷ Lu, and ¹³¹ I	xSPECT Quant ^{99m} Tc, ¹²³ I, ¹¹¹ In, ¹⁷⁷ Lu, and ¹³¹ I	
Quantitative SPECT studies for common SPECT radiopharmaceuticals in combination with all collimators	Broad Quantification (in units of Bq/ml, SUV, or count-rate-per-voxel)	Broad Quantification (in units of Bq/ml, SUV, or count-rate-per-voxel)	Broad Quantification (in units of Bq/ml, SUV, or count-rate-per-voxel)

xSPECT advanced specifications

xSPECT Quant: Accuracy of Bq/ml quantification ¹⁸ calibration method: NIST-traceable source		Symbia Pro.specta X7	Symbia Pro.specta X3	
Isotope/collimator	Uncertainty (95% confidence)	Uncertainty (95% confidence)		
^{99m} Tc LEHR	≤5%	≤5%		
^{99m} Tc LPHR	≤10%	≤10%		
¹²³ I LPHR	≤10%	≤10%		
¹²³ I MELP	≤10%	≤10%		
¹¹¹ In MELP	≤10%	≤10%		
¹³¹ I HE	≤10%	≤10%		
¹⁷⁷ Lu MELP	≤5%	≤5%		
¹⁷⁷ Lu MELP at 310 kcps incident count rate ¹⁰	≤10%	≤10%		
¹³¹ I at 310 kcps incident count rate ¹⁰	≤10%	≤10%		

Broad Quantification ¹⁹ : Reproducibility of Bq/ml quantification ²⁰ calibration method: Dose calibrator source measurement		Symbia Pro.specta X7	Symbia Pro.specta X3	Symbia Pro.specta Q3
Isotope/collimator	Reproducible within	Reproducible within	Reproducible within	Reproducible within
¹³¹ I HE	≤10%	≤10%	≤10%	≤10%
⁶⁷ Ga MELP	≤10%	≤10%	≤10%	≤10%

Room layout²¹



Scanner room size	3.74 m (12 ft 3.2 in) x 6.284 m (20 ft 7.4 in) recommended 3.74 m (12 ft 3.2 in) x 5.714 m (18 ft 9.0 in) minimum
Ceiling height	2,438 mm (8 ft) recommended 2,388 mm (7 ft 10 in) minimum

Installation and quality control specifications

Room diagram label	Item name	Weight	Heat output
	Symbia Pro.specta™		
SG	SPECT gantry	2,410 kg (5,313 lb)	<3,413 BTU/h, <1.0 kW ²²
Q3	Q3-CT gantry	1,282.3 kg (2,827 lb)	6.2 kW, 21,155 BTU/HR
X3	X3-CT gantry	1,300 kg (2,867 lb)	6.2 kW, 21,155 BTU/HR
X7	X7-CT gantry	1,300.8 kg (2,867.8 lb)	7.8 kW, 26,615 BTU/HR
T	Symbia Pro.specta imaging table	837.2 kg (1,845.8 lb)	
R	Symbia Pro.specta rear PHS	183 kg (403 lb)	
IEC	Siemens Healthineers power distributor Integrated Electronics Cabinet (IEC)	194 kg (428 lb)	0.15 kW, 512 BTU/HR

Collimator type	Weight
LEAP	22.6 kg (49.8 lb)
LEUHR	28.0 kg (61.8 lb)
LEHR	22.1 kg (48.7 lb)
MELP	63.5 kg (140.1 lb)
LPHR	33.6 kg (74 lb)
Pinhole	80.3 kg (177.0 lb)
HE	124.7 kg (275.0 lb)
SMARTZOOM	47.2 kg (104.0 lb)
SMARTZOOM HRX	56.1 kg (123.7 lb)
ICC housing	82 kg (182 lb)
ACC	45 kg (100 lb)
AQC	23 kg (50 lb)
Collimator cart	240 kg (528 lb)

Installation and quality control specifications

Power requirements	
Input voltage	Three-phase 380/400/420/440/460/480 VAC ~50/60 Hz
Electrical supply	For Symbia Pro.specta Q3 and X3 50 kVA ²²
	For Symbia Pro.specta X7 115 kVA ²²
Environment	
Floor loading	5.1 kg/sq cm (78 lb/sq in) maximum under the gantry
Ambient operating temperature	18-25°C (65-77°F)
Allowable temperature change	4.4°C (8°F) per hour
Humidity range	20-75% non-condensing
Maximum altitude	2,000 m a.s.l.
Standard quality control procedures	
Nuclear medicine	
Daily	Intrinsic verification or extrinsic verification
Weekly	Tuning
Monthly	– Intrinsic calibration and verification – Head alignment calibration with trending collimator
Every 6 months or per regulatory/license requirements (if applicable)	Sensitivity calibrations Leak test of the automated quality control device sources
Yearly	Head alignment calibration with every collimator
Computed tomography	
Daily	CT checkup at the start of the day, CT quality (DailyQA) at the start of the day, CT calibration 1 hour after checkup or if ring artifacts occur
Monthly	CT constancy test

Footnotes

- ¹ Available only for Symbia Pro.specta Q3
- ² Not available for Symbia Pro.specta Q3
- ³ Available only for Symbia Pro.specta X7
- ⁴ Distance from center of SPECT FOV to center of CT FOV
- ⁵ All collimators may not be supported by all detector configurations
- ⁶ For any point on the pallet at maximum 183 cm (6 ft) from the detector while the detector is at 25.4 cm (10 in) radial position
- ⁷ Distance from the edge of the detector housing to the edge of the FOV
- ⁸ Values are determined at the manufacturer's facility using methods described in NEMA Standards Publications NU 1-2018 "Performance measurements of Gamma Cameras."
- ⁹ With TrueCalc high count rate detector technology
- ¹⁰ Incident count rate
- ¹¹ Values measured in accordance with NEMA Standards Publication NU-1 2018 using 3/8" crystal
- ¹² Values measured using a 5-cm diameter phantom
- ¹³ Values measured with lines spaced 2 cm apart at the center of the collimator
- ¹⁴ Values measured in accordance with NEMA Standards Publication NU-1 2018 using 3/8" crystal. Sensitivities for pinhole collimators measured using a 9-cm diameter phantom. Resolution for pinhole collimator measured using a line separation of 6 cm (4-mm and 6-mm apertures) and 4 cm (8-mm aperture).
- ¹⁵ With HD FoV option. The image quality for the area outside the 50-cm scan field of view does not meet the image quality of the area inside the 50-cm scan field of view. Image artifacts may appear, depending on the patient setup and anatomy scanned. HD FoV cannot be used for scan FOV smaller than 50 cm.
- ¹⁶ Optional
- ¹⁷ PMMA phantom. Absorbed dose for reference material air. Maximum deviation $\pm 30\%$. Expected deviation $\pm 15\%$. Slice > 1 mm. Please note that these specifications are CTDI100 values.
- ¹⁸ Measured per NEMA NU1-2018 using a uniform cylinder phantom
- ¹⁹ Broad Quantification may be used for multiple SPECT isotopes beyond those listed.
- ²⁰ Using a uniform cylinder phantom
- ²¹ Example layout. Please request site-specific plans for your project.
- ²² Maximum power consumption during CT operation

Siemens Healthineers Molecular Imaging reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Sales representative for the most current information. Some options and functionality will not be available immediately on product release. Where certain options and functionality are not available on delivery, these will be delivered as part of subsequent software or hardware releases. Please confirm availability and timing with your representative.

This document provides information regarding technical specifications and standard and optional features. The listed specifications and features do not apply to all products. Note: technical data provided in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Symbia Pro.specta and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Trademarks and service marks used in this material are property of Siemens Medical Solutions USA or Siemens Healthcare GmbH. All other company, brand, product, and service names are the property of their respective holders.

“Siemens Healthineers” is considered a brand name. Its use is not intended to represent the legal entity to which this product is registered. Please contact your local Siemens Healthineers organization for further details.

All photographs © 2022 Siemens Healthcare GmbH.
All rights reserved.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Global Business Line

Siemens Medical Solutions USA, Inc.
Molecular Imaging
2501 North Barrington Road
Hoffman Estates, IL 60192
USA
Phone: +1 847 304-7700
siemens-healthineers.com/mi

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Molecular Imaging
2501 North Barrington Road
Hoffman Estates, IL 60192
USA
Phone: +1 847 304-7700
siemens-healthineers.com/mi