

MI-SPECT Software Versions

Release	Med Ver	e.soft ver	e.cam SNAC	Symbia SNAC	Brand IT	Syngo Version	Somaris Version
—	—	1.0	3.0	—	—	—	—
—	—	1.0	3.0.1	—	—	—	—
—	—	1.0	3.0.2	—	—	—	—
—	—	1.5	3.5	—	—	—	—
—	—	1.6	3.6	—	—	—	—
—	—	1.7	3.7	—	—	—	—
—	—	2.0	4.0	—	—	VB10D	—
—	—	2.1	4.1	—	—	VB10D	—
—	—	2.5	4.5	—	—	VD20H	—
—	—	3.0	5.0	—	—	VD20J	—
—	—	3.5	5.5	—	—	VD20K	—
—	—	4.0	6.0	—	2004B	VD20O	—
1.0	VA10	5.0	—	7.0	2005A	VD30A	VB20
1.1	VA11	5.1	—	7.1	2005B	VD30B	VB20
1.5	VA15	5.5	—	7.5	2005C	VD30A	VB20
2.0	VA20	6.0	8.0	8.0	2006A	VE31E (Nuc) VD30C (CT)	VB20
2.5	VA25	6.5	8.5	8.5	2007A	VE31F (Nuc) VD30C (CT)	VB20
3.0	VA30	7.0	—	9.0	2007B	VE31F	VB38
3.5	VA40	6.6 (e.cam)	8.6	—	2008A	VE31F	VB38
		7.5	—	9.5			
4.0	VA50A	8.0	—	10.0	2009A	VE32B (Nuc)	VB38
						VE31F (CT)	
4.1	VA50B	8.1	—	10.1	2009A	VE32B (Nuc)	VB38
						VE31F (CT)	
4.2	VA50C	8.2	—	10.2	2009A	VE32B (Nuc)	VB38
						VE31F (CT)	
4.5	VA60A	8.5	—	10.5	VA60A	VE32B (Nuc)	VB38
	VA46A	6.7 (e.cam)	8.7	—	VA46A	VE31F (CT)	
4.5	VA60B	8.5	—	10.5	VA60B	VE31F (e.cam)	VB38
	VA46B	6.7 (e.cam)	8.7	—	VA46B	VE32B (Nuc)	
4.5	VA60C	8.5	—	10.5	VA60C	VE31F (CT)	VB38
	VA46C	6.7 (e.cam)	8.7	—	VA46C	VE31F (e.cam)	
4.5	VA60E	8.5	—	12.1	VA60E	VE32B (Nuc)	VB38
	VA46D	6.7 (e.cam)	8.9	—	VA46D	VE31F (CT)	
4.5 SP2	VA61A	8.5	—	10.5	VA61A	VE31F (e.cam)	VB38
						VE32B (Nuc)	
Symbia T Series Only							
4.5 SP2	VA61B	8.5	—	10.5	VA61B	VE32B (Nuc)	VB38
						VE31F (CT)	
4.5 SP2	VA61C	8.5	—	10.5	VA61C	VE32B (Nuc)	VB38
						VE31F (CT)	
4.5 SP2	VA61D	8.5	—	12.1	VA61D	VE32B (Nuc)	VB38
						VE31F (CT)	
Symbia T Series Only							
Acq Cut	VA63A	8.6	—	10.6	VA63A	VE31F (CT)	VB38
5.0	VA70A	9.0	—	11.0	VA70A	VX70F (Nuc)	VB38
						VX68F (CT)	
						VX91B (DRS)	
6.0	VB10A	10.0	12.0	12.0	VB10A	VH10	VC20
6.0	VB10B	10.0	12.1	12.1	VB10B	VH10	VC20