

THEORY

Before an image can be transferred via DICOM, the sending computer must make a connection (an **association**) with the receiving computer, and then the two systems must agree (**negotiate**) the types of items that can be sent AND what syntax will be utilized for each type of item. As an example, a Siemens ESOFTE and a PACS first make an association, and then the two computers might negotiate that nuclear medicine studies will be sent using Implicit Little Endian syntax.

Little Endian means the least significant byte is sent first.

Big Endian means the most significant byte is sent first.

Implicit means the Value Representation of the attributes is not sent.

Explicit means the Value Representation of the attributes is sent.

In most DICOM applications, Implicit Little Endian is the default transfer syntax.

In DICOM terminology, this process is called **association/negotiation**.

To make this process easier to understand, here is an analogy:

Tom wants to send Gary a digital camera. First, Tom must get Gary's attention, i.e., make an association. Second, Tom tells Gary what he wants to send AND all the methods it can be sent. Then Gary responds with the shipping method he will accept. This second part is the negotiation.

Since Tom ships everything from software to electronics to jewelry, etc., Tom uses a standard introduction message detailing all the items he can ship, along with all their possible shipping methods. Then Gary replies with the shipping method he will accept for each type of item.

For example...

Tom: Hello Gary, my name is Tom.
I can ship software via...

- Internet
- FedEx
- UPS
- US Mail

I can ship electronics via...

- FedEx
- UPS

I can ship jewelry via...

- FedEx
- UPS

etc, etc, etc

Gary: Hello Tom, my name is Gary.
I will accept software via the internet.
I will accept electronics via UPS.
I will accept jewelry via FedEx.
etc, etc, etc

Gary might even refuse some offerings.

Tom: I can ship furniture via...

- BAX
- Allied

Gary: I will not accept furniture.

Once the association has been accomplished and the two parties have MUTUALLY agreed upon what types of items can shipped and MUTUALLY agreed upon a shipping method for each type of item, Tom announces exactly what he wants to send.

Tom: Gary, I am going to send you electronics via UPS.

Then Tom attempts to ship the digital camera to Gary via UPS.

DICOM works basically the same way.

In the following example, a Siemens ESOFTE is attempting to transfer a nuclear medicine study to a PACS. To begin, the ESOFTE (**AET=ESOFTE**) has called the PACS (**AET=AQNET_AE**) on port 105. ESOFTE then announced the 74 different types of studies it could transfer (e.g., CT, MR, NM, US, etc.) and listed all the available transfer syntaxes for each. For example, for nuclear medicine studies, ESOFTE offered...

- JPEG Lossless
- JPEG Baseline
- Implicit Little Endian
- Explicit Little Endian
- Explicit Big Endian

A Siemens service engineer will probably recognize these five possible transfer syntaxes from the Options → Service → Local Service → Configuration page with *syngo*. In the Local Service → Configuration interface, a service engineer has the ability of enabling or disabling any of these syntaxes on a node by node basis.

In the following example, AQNET_AE responded with the transfer syntax it would accept for each type of study. **Note that AQNET_AE rejected many of the proposed transfer syntaxes. Many people mistakenly interpret these rejections as a problem; however, as long as the two systems can agree on a syntax for the item being transferred, the negotiation is successful.** Using the previous digital camera analogy, it doesn't matter that Gary rejected Tom's offer to ship furniture via BAX or Allied, because Tom was not sending furniture.

In the following example, AQNET_AE indicated it would only accept Implicit Little Endian for nuclear medicine studies; therefore, ESOFTE should only attempt to send a NM study using that syntax.

Next, ESOFTE stated that it wanted AQNET_AE to store a standard NM study. Furthermore, ESOFTE indicated the study would be sent using Implicit Little Endian syntax. Afterwards, the file transfer commenced. Next, the study was received by AQNET_AE. The association/negotiation was successful - and the study transfer was also a success!

The following capture log is an excellent example because it allows someone to view a DICOM transfer when everything happened exactly as planned. Comments have been added to make reading the log easier.

START WITH THE SENDER CALLING THE RECEIVING AET

```
(5528) 06-26 13:20:46.09 MC3 T3: REQUEST Association Parameters:
(5528) 06-26 13:20:46.09 MC3 T3: | Calling_presentation_address: ESOFTB
(5528) 06-26 13:20:46.09 MC3 T3: | Called_presentation_address: aqnet_ae:105
(5528) 06-26 13:20:46.09 MC3 T3: | Called_application_title: AQNET_AE
(5528) 06-26 13:20:46.09 MC3 T3: | Calling_application_title: ESOFTB
```

HERE IS THE STANDARD ACN

```
(5528) 06-26 13:20:46.09 MC3 T3: | Application_context_name:
1.2.840.10008.3.1.1.1
```

MAX PROTOCOL DATA UNIT FOR THE SENDER (ESOFT) IS 28672

```
(5528) 06-26 13:20:46.10 MC3 T3: | PDU_maximum_length: 28672
```

THIS IS AN SIEMENS ESOFT IMPLEMENTATION

```
(5528) 06-26 13:20:46.10 MC3 T3: | Implem_class_uid:
1.3.12.2.1107.5.9.20000101
```

VD20K IS ESOFT 3.5

```
(5528) 06-26 13:20:46.10 MC3 T3: | Version_name: SIEMENS_SWFVD20K
(5528) 06-26 13:20:46.10 MC3 T3: | Max_operations_invoked: Not Sent
(5528) 06-26 13:20:46.10 MC3 T3: | Max_operations_performed: Not Sent
```

THE SENDER TELLS THE RECEIVER WHICH "TRANSFER SYNTAXES" THE SENDER CAN OFFER.

NM IS PRES CONTEXT #13 & #15.

```
(5528) 06-26 13:20:46.10 MC3 T3: | PRESENTATION CONTEXTS:
(5528) 06-26 13:20:46.10 MC3 T3: | (1)Pres_context_id: 1
(5528) 06-26 13:20:46.10 MC3 T3: | (1)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1 (STANDARD_CR)
(5528) 06-26 13:20:46.10 MC3 T3: | (1)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
```

(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | (1)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.10 MC3 T3: | (1)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.10 MC3 T3: | (1)No extended application info
(5528) 06-26 13:20:46.10 MC3 T3: | (2)Pres_context_id: 3
(5528) 06-26 13:20:46.10 MC3 T3: | (2)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1 (STANDARD_CR)
(5528) 06-26 13:20:46.10 MC3 T3: | (2)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | (2)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.10 MC3 T3: | (2)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.10 MC3 T3: | (2)No extended application info
(5528) 06-26 13:20:46.10 MC3 T3: | (3)Pres_context_id: 5
(5528) 06-26 13:20:46.10 MC3 T3: | (3)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.2 (STANDARD_CT)
(5528) 06-26 13:20:46.10 MC3 T3: | (3)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | (3)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.10 MC3 T3: | (3)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.10 MC3 T3: | (3)No extended application info
(5528) 06-26 13:20:46.10 MC3 T3: | (4)Pres_context_id: 7
(5528) 06-26 13:20:46.10 MC3 T3: | (4)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.2 (STANDARD_CT)
(5528) 06-26 13:20:46.10 MC3 T3: | (4)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little

Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | (4)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.10 MC3 T3: | (4)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.10 MC3 T3: | (4)No extended application info
(5528) 06-26 13:20:46.10 MC3 T3: | (5)Pres_context_id: 9
(5528) 06-26 13:20:46.10 MC3 T3: | (5)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.4 (STANDARD_MR)
(5528) 06-26 13:20:46.10 MC3 T3: | (5)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.10 MC3 T3: | (5)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.10 MC3 T3: | (5)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.10 MC3 T3: | (5)No extended application info
(5528) 06-26 13:20:46.10 MC3 T3: | (6)Pres_context_id: 11
(5528) 06-26 13:20:46.10 MC3 T3: | (6)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.4 (STANDARD_MR)
(5528) 06-26 13:20:46.10 MC3 T3: | (6)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.10 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | (6)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.12 MC3 T3: | (6)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.12 MC3 T3: | (6)No extended application info

**FOR NM, THE SENDER OFFERS...
JPEG LOSSLESS, JPEG BASELESS,
IMPLICIT LITTLE ENDIAN,**

EXPLICIT LITTLE ENDIAN, & EXPLICIT BIG ENDIAN.

```
(5528) 06-26 13:20:46.12 MC3 T3: | (7)Pres_context_id: 13
(5528) 06-26 13:20:46.12 MC3 T3: | (7)Abstract_syntax:
(5528) 06-26 13:20:46.12 MC3 T3: | (7)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.20 (STANDARD_NM)
(5528) 06-26 13:20:46.12 MC3 T3: | (7)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | (7)SCU Role: Not sent (requestor supports
      SCU)
(5528) 06-26 13:20:46.12 MC3 T3: | (7)SCP Role: Not sent (requestor supports
      SCP)
(5528) 06-26 13:20:46.12 MC3 T3: | (7)No extended application info
(5528) 06-26 13:20:46.12 MC3 T3: | (8)Pres_context_id: 15
(5528) 06-26 13:20:46.12 MC3 T3: | (8)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.20 (STANDARD_NM)
(5528) 06-26 13:20:46.12 MC3 T3: | (8)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
      4)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | (8)SCU Role: Not sent (requestor supports
      SCU)
(5528) 06-26 13:20:46.12 MC3 T3: | (8)SCP Role: Not sent (requestor supports
      SCP)
(5528) 06-26 13:20:46.12 MC3 T3: | (8)No extended application info
(5528) 06-26 13:20:46.12 MC3 T3: | (9)Pres_context_id: 17
(5528) 06-26 13:20:46.12 MC3 T3: | (9)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.7 (STANDARD_SEC_CAPTURE)
(5528) 06-26 13:20:46.12 MC3 T3: | (9)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
      Endian)
```

(5528) 06-26 13:20:46.12 MC3 T3: | (9)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.12 MC3 T3: | (9)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.12 MC3 T3: | (9)No extended application info
(5528) 06-26 13:20:46.12 MC3 T3: | (10)Pres_context_id: 19
(5528) 06-26 13:20:46.12 MC3 T3: | (10)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.7 (STANDARD_SEC_CAPTURE)
(5528) 06-26 13:20:46.12 MC3 T3: | (10)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | (10)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.12 MC3 T3: | (10)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.12 MC3 T3: | (10)No extended application info
(5528) 06-26 13:20:46.12 MC3 T3: | (11)Pres_context_id: 21
(5528) 06-26 13:20:46.12 MC3 T3: | (11)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.6.1 (STANDARD_US)
(5528) 06-26 13:20:46.12 MC3 T3: | (11)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | (11)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.12 MC3 T3: | (11)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.12 MC3 T3: | (11)No extended application info
(5528) 06-26 13:20:46.12 MC3 T3: | (12)Pres_context_id: 23
(5528) 06-26 13:20:46.12 MC3 T3: | (12)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.6.1 (STANDARD_US)
(5528) 06-26 13:20:46.12 MC3 T3: | (12)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.12 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | (12)SCU Role: Not sent (requestor

supports SCU)
(5528) 06-26 13:20:46.14 MC3 T3: | (12)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.14 MC3 T3: | (12)No extended application info
(5528) 06-26 13:20:46.14 MC3 T3: | (13)Pres_context_id: 25
(5528) 06-26 13:20:46.14 MC3 T3: | (13)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.3.1 (STANDARD_US_MF)
(5528) 06-26 13:20:46.14 MC3 T3: | (13)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | (13)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.14 MC3 T3: | (13)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.14 MC3 T3: | (13)No extended application info
(5528) 06-26 13:20:46.14 MC3 T3: | (14)Pres_context_id: 27
(5528) 06-26 13:20:46.14 MC3 T3: | (14)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.3.1 (STANDARD_US_MF)
(5528) 06-26 13:20:46.14 MC3 T3: | (14)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | (14)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.14 MC3 T3: | (14)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.14 MC3 T3: | (14)No extended application info
(5528) 06-26 13:20:46.14 MC3 T3: | (15)Pres_context_id: 29
(5528) 06-26 13:20:46.14 MC3 T3: | (15)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.128 (STANDARD_PET)
(5528) 06-26 13:20:46.14 MC3 T3: | (15)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | (15)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.14 MC3 T3: | (15)SCP Role: Not sent (requestor


```
supports SCP)
(5528) 06-26 13:20:46.14 MC3 T3: | (15)No extended application info
(5528) 06-26 13:20:46.14 MC3 T3: | (16)Pres_context_id: 31
(5528) 06-26 13:20:46.14 MC3 T3: | (16)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.128 (STANDARD_PET)
(5528) 06-26 13:20:46.14 MC3 T3: | (16)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.4.51 (JPEG Extended 2
      4)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | (16)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.14 MC3 T3: | (16)SCP Role: Not sent (requestor
      supports SCP)
(5528) 06-26 13:20:46.14 MC3 T3: | (16)No extended application info
(5528) 06-26 13:20:46.14 MC3 T3: | (17)Pres_context_id: 33
(5528) 06-26 13:20:46.14 MC3 T3: | (17)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.1.1 (STANDARD_WAVEFORM_12_LEAD_ECG)
(5528) 06-26 13:20:46.14 MC3 T3: | (17)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: | (17)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.14 MC3 T3: | (17)SCP Role: Not sent (requestor
      supports SCP)
(5528) 06-26 13:20:46.14 MC3 T3: | (17)No extended application info
(5528) 06-26 13:20:46.14 MC3 T3: | (18)Pres_context_id: 35
(5528) 06-26 13:20:46.14 MC3 T3: | (18)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.1.1 (STANDARD_WAVEFORM_12_LEAD_ECG)
(5528) 06-26 13:20:46.14 MC3 T3: | (18)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.4.51 (JPEG Extended 2
      4)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.14 MC3 T3: |   1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | (18)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.15 MC3 T3: | (18)SCP Role: Not sent (requestor
      supports SCP)
```

(5528) 06-26 13:20:46.15 MC3 T3: | (18)No extended application info
(5528) 06-26 13:20:46.15 MC3 T3: | (19)Pres_context_id: 37
(5528) 06-26 13:20:46.15 MC3 T3: | (19)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.2 (STANDARD_WAVEFORM_GENERAL_ECG)
(5528) 06-26 13:20:46.15 MC3 T3: | (19)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | (19)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.15 MC3 T3: | (19)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.15 MC3 T3: | (19)No extended application info
(5528) 06-26 13:20:46.15 MC3 T3: | (20)Pres_context_id: 39
(5528) 06-26 13:20:46.15 MC3 T3: | (20)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.2 (STANDARD_WAVEFORM_GENERAL_ECG)
(5528) 06-26 13:20:46.15 MC3 T3: | (20)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | (20)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.15 MC3 T3: | (20)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.15 MC3 T3: | (20)No extended application info
(5528) 06-26 13:20:46.15 MC3 T3: | (21)Pres_context_id: 41
(5528) 06-26 13:20:46.15 MC3 T3: | (21)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.3
(STANDARD_WAVEFORM_AMBULATORY_ECG)
(5528) 06-26 13:20:46.15 MC3 T3: | (21)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | (21)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.15 MC3 T3: | (21)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.15 MC3 T3: | (21)No extended application info

(5528) 06-26 13:20:46.15 MC3 T3: | (22)Pres_context_id: 43
(5528) 06-26 13:20:46.15 MC3 T3: | (22)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.3
(STANDARD_WAVEFORM_AMBULATORY_ECG)
(5528) 06-26 13:20:46.15 MC3 T3: | (22)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | (22)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.15 MC3 T3: | (22)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.15 MC3 T3: | (22)No extended application info
(5528) 06-26 13:20:46.15 MC3 T3: | (23)Pres_context_id: 45
(5528) 06-26 13:20:46.15 MC3 T3: | (23)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.2.1 (STANDARD_WAVEFORM_HEMODYNAMIC)
(5528) 06-26 13:20:46.15 MC3 T3: | (23)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.15 MC3 T3: | (23)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.15 MC3 T3: | (23)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.15 MC3 T3: | (23)No extended application info
(5528) 06-26 13:20:46.15 MC3 T3: | (24)Pres_context_id: 47
(5528) 06-26 13:20:46.15 MC3 T3: | (24)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.2.1 (STANDARD_WAVEFORM_HEMODYNAMIC)
(5528) 06-26 13:20:46.15 MC3 T3: | (24)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.15 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | (24)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.17 MC3 T3: | (24)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.17 MC3 T3: | (24)No extended application info

(5528) 06-26 13:20:46.17 MC3 T3: | (25)Pres_context_id: 49
(5528) 06-26 13:20:46.17 MC3 T3: | (25)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.3.1 (STANDARD_WAVEFORM_CARDIAC_EP)
(5528) 06-26 13:20:46.17 MC3 T3: | (25)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | (25)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.17 MC3 T3: | (25)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.17 MC3 T3: | (25)No extended application info
(5528) 06-26 13:20:46.17 MC3 T3: | (26)Pres_context_id: 51
(5528) 06-26 13:20:46.17 MC3 T3: | (26)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.3.1 (STANDARD_WAVEFORM_CARDIAC_EP)
(5528) 06-26 13:20:46.17 MC3 T3: | (26)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | (26)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.17 MC3 T3: | (26)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.17 MC3 T3: | (26)No extended application info
(5528) 06-26 13:20:46.17 MC3 T3: | (27)Pres_context_id: 53
(5528) 06-26 13:20:46.17 MC3 T3: | (27)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.4.1
(STANDARD_WAVEFORM_BASIC_VOICE_AU)
(5528) 06-26 13:20:46.17 MC3 T3: | (27)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | (27)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.17 MC3 T3: | (27)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.17 MC3 T3: | (27)No extended application info
(5528) 06-26 13:20:46.17 MC3 T3: | (28)Pres_context_id: 55

(5528) 06-26 13:20:46.17 MC3 T3: | (28)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.4.1
(STANDARD_WAVEFORM_BASIC_VOICE_AU)

(5528) 06-26 13:20:46.17 MC3 T3: | (28)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.17 MC3 T3: | (28)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.17 MC3 T3: | (28)SCP Role: Not sent (requestor
supports SCP)

(5528) 06-26 13:20:46.17 MC3 T3: | (28)No extended application info
(5528) 06-26 13:20:46.17 MC3 T3: | (29)Pres_context_id: 57
(5528) 06-26 13:20:46.17 MC3 T3: | (29)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.1 (STANDARD_XRAY_ANGIO)

(5528) 06-26 13:20:46.17 MC3 T3: | (29)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.17 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.17 MC3 T3: | (29)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.17 MC3 T3: | (29)SCP Role: Not sent (requestor
supports SCP)

(5528) 06-26 13:20:46.18 MC3 T3: | (29)No extended application info
(5528) 06-26 13:20:46.18 MC3 T3: | (30)Pres_context_id: 59
(5528) 06-26 13:20:46.18 MC3 T3: | (30)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.1 (STANDARD_XRAY_ANGIO)

(5528) 06-26 13:20:46.18 MC3 T3: | (30)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.18 MC3 T3: | (30)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.18 MC3 T3: | (30)SCP Role: Not sent (requestor
supports SCP)

(5528) 06-26 13:20:46.18 MC3 T3: | (30)No extended application info
(5528) 06-26 13:20:46.18 MC3 T3: | (31)Pres_context_id: 61

(5528) 06-26 13:20:46.18 MC3 T3: | (31)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.3 (STANDARD_XRAY_ANGIO_BIPLANE)
(5528) 06-26 13:20:46.18 MC3 T3: | (31)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | (31)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.18 MC3 T3: | (31)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.18 MC3 T3: | (31)No extended application info
(5528) 06-26 13:20:46.18 MC3 T3: | (32)Pres_context_id: 63
(5528) 06-26 13:20:46.18 MC3 T3: | (32)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.3 (STANDARD_XRAY_ANGIO_BIPLANE)
(5528) 06-26 13:20:46.18 MC3 T3: | (32)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | (32)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.18 MC3 T3: | (32)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.18 MC3 T3: | (32)No extended application info
(5528) 06-26 13:20:46.18 MC3 T3: | (33)Pres_context_id: 65
(5528) 06-26 13:20:46.18 MC3 T3: | (33)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.2 (STANDARD_XRAY_RF)
(5528) 06-26 13:20:46.18 MC3 T3: | (33)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | (33)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.18 MC3 T3: | (33)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.18 MC3 T3: | (33)No extended application info
(5528) 06-26 13:20:46.18 MC3 T3: | (34)Pres_context_id: 67
(5528) 06-26 13:20:46.18 MC3 T3: | (34)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.2 (STANDARD_XRAY_RF)

(5528) 06-26 13:20:46.18 MC3 T3: | (34)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.18 MC3 T3: | (34)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.18 MC3 T3: | (34)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.18 MC3 T3: | (34)No extended application info
(5528) 06-26 13:20:46.18 MC3 T3: | (35)Pres_context_id: 69
(5528) 06-26 13:20:46.18 MC3 T3: | (35)Abstract_syntax: 1.3.12.2.1107.5.9.1
(PRIVATE_MEDCOM_CSA_NON_IMAGE)
(5528) 06-26 13:20:46.20 MC3 T3: | (35)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.20 MC3 T3: | (35)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.20 MC3 T3: | (35)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.20 MC3 T3: | (35)No extended application info
(5528) 06-26 13:20:46.20 MC3 T3: | (36)Pres_context_id: 71
(5528) 06-26 13:20:46.20 MC3 T3: | (36)Abstract_syntax: 1.3.12.2.1107.5.9.1
(PRIVATE_MEDCOM_CSA_NON_IMAGE)
(5528) 06-26 13:20:46.20 MC3 T3: | (36)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.20 MC3 T3: | (36)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.20 MC3 T3: | (36)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.20 MC3 T3: | (36)No extended application info
(5528) 06-26 13:20:46.20 MC3 T3: | (37)Pres_context_id: 73
(5528) 06-26 13:20:46.20 MC3 T3: | (37)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.3 (STANDARD_RT_STRUCTURE_SET)
(5528) 06-26 13:20:46.20 MC3 T3: | (37)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | (37)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.20 MC3 T3: | (37)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.20 MC3 T3: | (37)No extended application info

(5528) 06-26 13:20:46.20 MC3 T3: | (38)Pres_context_id: 75

(5528) 06-26 13:20:46.20 MC3 T3: | (38)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.3 (STANDARD_RT_STRUCTURE_SET)

(5528) 06-26 13:20:46.20 MC3 T3: | (38)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | (38)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.20 MC3 T3: | (38)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.20 MC3 T3: | (38)No extended application info

(5528) 06-26 13:20:46.20 MC3 T3: | (39)Pres_context_id: 77

(5528) 06-26 13:20:46.20 MC3 T3: | (39)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.1 (STANDARD_DX_PRESENT)

(5528) 06-26 13:20:46.20 MC3 T3: | (39)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.20 MC3 T3: | (39)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.20 MC3 T3: | (39)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.20 MC3 T3: | (39)No extended application info

(5528) 06-26 13:20:46.20 MC3 T3: | (40)Pres_context_id: 79

(5528) 06-26 13:20:46.20 MC3 T3: | (40)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.1 (STANDARD_DX_PRESENT)

(5528) 06-26 13:20:46.20 MC3 T3: | (40)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.20 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | (40)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.21 MC3 T3: | (40)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.21 MC3 T3: | (40)No extended application info

(5528) 06-26 13:20:46.21 MC3 T3: | (41)Pres_context_id: 81

(5528) 06-26 13:20:46.21 MC3 T3: | (41)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.1.1 (STANDARD_DX_PROCESS)

(5528) 06-26 13:20:46.21 MC3 T3: | (41)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | (41)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.21 MC3 T3: | (41)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.21 MC3 T3: | (41)No extended application info

(5528) 06-26 13:20:46.21 MC3 T3: | (42)Pres_context_id: 83

(5528) 06-26 13:20:46.21 MC3 T3: | (42)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.1.1 (STANDARD_DX_PROCESS)

(5528) 06-26 13:20:46.21 MC3 T3: | (42)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.21 MC3 T3: | (42)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.21 MC3 T3: | (42)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.21 MC3 T3: | (42)No extended application info

(5528) 06-26 13:20:46.21 MC3 T3: | (43)Pres_context_id: 85

(5528) 06-26 13:20:46.21 MC3 T3: | (43)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.3 (STANDARD_IO_PRESENT)

(5528) 06-26 13:20:46.21 MC3 T3: | (43)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)

(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.21 MC3 T3: | (43)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.21 MC3 T3: | (43)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.21 MC3 T3: | (43)No extended application info
(5528) 06-26 13:20:46.21 MC3 T3: | (44)Pres_context_id: 87
(5528) 06-26 13:20:46.21 MC3 T3: | (44)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.3 (STANDARD_IO_PRESENT)
(5528) 06-26 13:20:46.21 MC3 T3: | (44)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.21 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.21 MC3 T3: | (44)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.21 MC3 T3: | (44)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.21 MC3 T3: | (44)No extended application info
(5528) 06-26 13:20:46.21 MC3 T3: | (45)Pres_context_id: 89
(5528) 06-26 13:20:46.21 MC3 T3: | (45)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.3.1 (STANDARD_IO_PROCESS)
(5528) 06-26 13:20:46.21 MC3 T3: | (45)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | (45)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.23 MC3 T3: | (45)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.23 MC3 T3: | (45)No extended application info
(5528) 06-26 13:20:46.23 MC3 T3: | (46)Pres_context_id: 91
(5528) 06-26 13:20:46.23 MC3 T3: | (46)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.3.1 (STANDARD_IO_PROCESS)
(5528) 06-26 13:20:46.23 MC3 T3: | (46)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little

Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | (46)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.23 MC3 T3: | (46)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.23 MC3 T3: | (46)No extended application info
(5528) 06-26 13:20:46.23 MC3 T3: | (47)Pres_context_id: 93
(5528) 06-26 13:20:46.23 MC3 T3: | (47)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2 (STANDARD_MG_PRESENT)
(5528) 06-26 13:20:46.23 MC3 T3: | (47)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | (47)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.23 MC3 T3: | (47)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.23 MC3 T3: | (47)No extended application info
(5528) 06-26 13:20:46.23 MC3 T3: | (48)Pres_context_id: 95
(5528) 06-26 13:20:46.23 MC3 T3: | (48)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2 (STANDARD_MG_PRESENT)
(5528) 06-26 13:20:46.23 MC3 T3: | (48)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | (48)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.23 MC3 T3: | (48)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.23 MC3 T3: | (48)No extended application info
(5528) 06-26 13:20:46.23 MC3 T3: | (49)Pres_context_id: 97
(5528) 06-26 13:20:46.23 MC3 T3: | (49)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2.1 (STANDARD_MG_PROCESS)
(5528) 06-26 13:20:46.23 MC3 T3: | (49)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little

Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.23 MC3 T3: | (49)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.23 MC3 T3: | (49)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.23 MC3 T3: | (49)No extended application info
(5528) 06-26 13:20:46.23 MC3 T3: | (50)Pres_context_id: 99
(5528) 06-26 13:20:46.23 MC3 T3: | (50)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2.1 (STANDARD_MG_PROCESS)
(5528) 06-26 13:20:46.25 MC3 T3: | (50)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | (50)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.25 MC3 T3: | (50)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.25 MC3 T3: | (50)No extended application info
(5528) 06-26 13:20:46.25 MC3 T3: | (51)Pres_context_id: 101
(5528) 06-26 13:20:46.25 MC3 T3: | (51)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.481.2 (STANDARD_RT_DOSE)
(5528) 06-26 13:20:46.25 MC3 T3: | (51)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | (51)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.25 MC3 T3: | (51)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.25 MC3 T3: | (51)No extended application info
(5528) 06-26 13:20:46.25 MC3 T3: | (52)Pres_context_id: 103
(5528) 06-26 13:20:46.25 MC3 T3: | (52)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.481.2 (STANDARD_RT_DOSE)
(5528) 06-26 13:20:46.25 MC3 T3: | (52)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)

(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | (52)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.25 MC3 T3: | (52)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.25 MC3 T3: | (52)No extended application info
(5528) 06-26 13:20:46.25 MC3 T3: | (53)Pres_context_id: 105
(5528) 06-26 13:20:46.25 MC3 T3: | (53)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.1 (STANDARD_RT_IMAGE)
(5528) 06-26 13:20:46.25 MC3 T3: | (53)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | (53)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.25 MC3 T3: | (53)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.25 MC3 T3: | (53)No extended application info
(5528) 06-26 13:20:46.25 MC3 T3: | (54)Pres_context_id: 107
(5528) 06-26 13:20:46.25 MC3 T3: | (54)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.1 (STANDARD_RT_IMAGE)
(5528) 06-26 13:20:46.25 MC3 T3: | (54)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2 4)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)
(5528) 06-26 13:20:46.25 MC3 T3: | (54)SCU Role: Not sent (requestor supports SCU)
(5528) 06-26 13:20:46.26 MC3 T3: | (54)SCP Role: Not sent (requestor supports SCP)
(5528) 06-26 13:20:46.26 MC3 T3: | (54)No extended application info
(5528) 06-26 13:20:46.26 MC3 T3: | (55)Pres_context_id: 109
(5528) 06-26 13:20:46.26 MC3 T3: | (55)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.5 (STANDARD_RT_PLAN)
(5528) 06-26 13:20:46.26 MC3 T3: | (55)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little Endian)
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | (55)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.26 MC3 T3: | (55)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.26 MC3 T3: | (55)No extended application info

(5528) 06-26 13:20:46.26 MC3 T3: | (56)Pres_context_id: 111

(5528) 06-26 13:20:46.26 MC3 T3: | (56)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.5 (STANDARD_RT_PLAN)

(5528) 06-26 13:20:46.26 MC3 T3: | (56)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | (56)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.26 MC3 T3: | (56)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.26 MC3 T3: | (56)No extended application info

(5528) 06-26 13:20:46.26 MC3 T3: | (57)Pres_context_id: 113

(5528) 06-26 13:20:46.26 MC3 T3: | (57)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.1.1 (PATIENT_ROOT_QR_FIND)

(5528) 06-26 13:20:46.26 MC3 T3: | (57)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | (57)SCU Role: Not sent (requestor supports SCU)

(5528) 06-26 13:20:46.26 MC3 T3: | (57)SCP Role: Not sent (requestor supports SCP)

(5528) 06-26 13:20:46.26 MC3 T3: | (57)No extended application info

(5528) 06-26 13:20:46.26 MC3 T3: | (58)Pres_context_id: 115

(5528) 06-26 13:20:46.26 MC3 T3: | (58)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.1.1 (PATIENT_ROOT_QR_FIND)

(5528) 06-26 13:20:46.26 MC3 T3: | (58)Proposed Transfer Syntaxes:

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.26 MC3 T3: | (58)SCU Role: Not sent (requestor supports SCU)

```
supports SCU)
(5528) 06-26 13:20:46.26 MC3 T3: | (58)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.26 MC3 T3: | (58)No extended application info
(5528) 06-26 13:20:46.26 MC3 T3: | (59)Pres_context_id: 117
(5528) 06-26 13:20:46.26 MC3 T3: | (59)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.1 (STUDY_ROOT_QR_FIND)
(5528) 06-26 13:20:46.26 MC3 T3: | (59)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.26 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | (59)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.28 MC3 T3: | (59)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.28 MC3 T3: | (59)No extended application info
(5528) 06-26 13:20:46.28 MC3 T3: | (60)Pres_context_id: 119
(5528) 06-26 13:20:46.28 MC3 T3: | (60)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.1 (STUDY_ROOT_QR_FIND)
(5528) 06-26 13:20:46.28 MC3 T3: | (60)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | (60)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.28 MC3 T3: | (60)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.28 MC3 T3: | (60)No extended application info
(5528) 06-26 13:20:46.28 MC3 T3: | (61)Pres_context_id: 121
(5528) 06-26 13:20:46.28 MC3 T3: | (61)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.3.1 (PATIENT_STUDY_ONLY_QR_FIND)
(5528) 06-26 13:20:46.28 MC3 T3: | (61)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | (61)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.28 MC3 T3: | (61)SCP Role: Not sent (requestor
```

```
supports SCP)
(5528) 06-26 13:20:46.28 MC3 T3: | (61)No extended application info
(5528) 06-26 13:20:46.28 MC3 T3: | (62)Pres_context_id: 123
(5528) 06-26 13:20:46.28 MC3 T3: | (62)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.3.1 (PATIENT_STUDY_ONLY_QR_FIND)
(5528) 06-26 13:20:46.28 MC3 T3: | (62)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.4.51 (JPEG Extended 2
      4)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | (62)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.28 MC3 T3: | (62)SCP Role: Not sent (requestor
      supports SCP)
(5528) 06-26 13:20:46.28 MC3 T3: | (62)No extended application info
(5528) 06-26 13:20:46.28 MC3 T3: | (63)Pres_context_id: 125
(5528) 06-26 13:20:46.28 MC3 T3: | (63)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.3 (PATIENT_ROOT_QR_GET)
(5528) 06-26 13:20:46.28 MC3 T3: | (63)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.28 MC3 T3: |   1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.28 MC3 T3: | (63)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.28 MC3 T3: | (63)SCP Role: Not sent (requestor
      supports SCP)
(5528) 06-26 13:20:46.28 MC3 T3: | (63)No extended application info
(5528) 06-26 13:20:46.29 MC3 T3: | (64)Pres_context_id: 127
(5528) 06-26 13:20:46.29 MC3 T3: | (64)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.3 (PATIENT_ROOT_QR_GET)
(5528) 06-26 13:20:46.29 MC3 T3: | (64)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.29 MC3 T3: |   1.2.840.10008.1.2.4.51 (JPEG Extended 2
      4)
(5528) 06-26 13:20:46.29 MC3 T3: |   1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.29 MC3 T3: |   1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.29 MC3 T3: |   1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.29 MC3 T3: |   1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | (64)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.29 MC3 T3: | (64)SCP Role: Not sent (requestor
      supports SCP)
```


(5528) 06-26 13:20:46.29 MC3 T3: | (64)No extended application info
(5528) 06-26 13:20:46.29 MC3 T3: | (65)Pres_context_id: 129
(5528) 06-26 13:20:46.29 MC3 T3: | (65)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.3.3 (PATIENT_STUDY_ONLY_QR_GET)
(5528) 06-26 13:20:46.29 MC3 T3: | (65)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | (65)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.29 MC3 T3: | (65)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.29 MC3 T3: | (65)No extended application info
(5528) 06-26 13:20:46.29 MC3 T3: | (66)Pres_context_id: 131
(5528) 06-26 13:20:46.29 MC3 T3: | (66)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.3.3 (PATIENT_STUDY_ONLY_QR_GET)
(5528) 06-26 13:20:46.29 MC3 T3: | (66)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | (66)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.29 MC3 T3: | (66)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.29 MC3 T3: | (66)No extended application info
(5528) 06-26 13:20:46.29 MC3 T3: | (67)Pres_context_id: 133
(5528) 06-26 13:20:46.29 MC3 T3: | (67)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.3 (STUDY_ROOT_QR_GET)
(5528) 06-26 13:20:46.29 MC3 T3: | (67)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.29 MC3 T3: | (67)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.29 MC3 T3: | (67)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.29 MC3 T3: | (67)No extended application info
(5528) 06-26 13:20:46.29 MC3 T3: | (68)Pres_context_id: 135

(5528) 06-26 13:20:46.29 MC3 T3: | (68)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.3 (STUDY_ROOT_QR_GET)

(5528) 06-26 13:20:46.29 MC3 T3: | (68)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)

(5528) 06-26 13:20:46.29 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | (68)SCU Role: Not sent (requestor
supports SCU)

(5528) 06-26 13:20:46.31 MC3 T3: | (68)SCP Role: Not sent (requestor
supports SCP)

(5528) 06-26 13:20:46.31 MC3 T3: | (68)No extended application info
(5528) 06-26 13:20:46.31 MC3 T3: | (69)Pres_context_id: 137

(5528) 06-26 13:20:46.31 MC3 T3: | (69)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.1.2 (PATIENT_ROOT_QR_MOVE)

(5528) 06-26 13:20:46.31 MC3 T3: | (69)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | (69)SCU Role: Not sent (requestor
supports SCU)

(5528) 06-26 13:20:46.31 MC3 T3: | (69)SCP Role: Not sent (requestor
supports SCP)

(5528) 06-26 13:20:46.31 MC3 T3: | (69)No extended application info
(5528) 06-26 13:20:46.31 MC3 T3: | (70)Pres_context_id: 139

(5528) 06-26 13:20:46.31 MC3 T3: | (70)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.1.2 (PATIENT_ROOT_QR_MOVE)

(5528) 06-26 13:20:46.31 MC3 T3: | (70)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)

(5528) 06-26 13:20:46.31 MC3 T3: | (70)SCU Role: Not sent (requestor
supports SCU)

(5528) 06-26 13:20:46.31 MC3 T3: | (70)SCP Role: Not sent (requestor
supports SCP)

(5528) 06-26 13:20:46.31 MC3 T3: | (70)No extended application info
(5528) 06-26 13:20:46.31 MC3 T3: | (71)Pres_context_id: 141

(5528) 06-26 13:20:46.31 MC3 T3: | (71)Abstract_syntax:

```
1.2.840.10008.5.1.4.1.2.3.2 (PATIENT_STUDY_ONLY_QR_MOVE)
(5528) 06-26 13:20:46.31 MC3 T3: | (71)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.31 MC3 T3: | (71)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.31 MC3 T3: | (71)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.31 MC3 T3: | (71)No extended application info
(5528) 06-26 13:20:46.31 MC3 T3: | (72)Pres_context_id: 143
(5528) 06-26 13:20:46.31 MC3 T3: | (72)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.3.2 (PATIENT_STUDY_ONLY_QR_MOVE)
(5528) 06-26 13:20:46.31 MC3 T3: | (72)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
4)
(5528) 06-26 13:20:46.31 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | (72)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.32 MC3 T3: | (72)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.32 MC3 T3: | (72)No extended application info
(5528) 06-26 13:20:46.32 MC3 T3: | (73)Pres_context_id: 145
(5528) 06-26 13:20:46.32 MC3 T3: | (73)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.2 (STUDY_ROOT_QR_MOVE)
(5528) 06-26 13:20:46.32 MC3 T3: | (73)Proposed Transfer Syntaxes:
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
Hier 14)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.4.50 (JPEG Baseline)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | (73)SCU Role: Not sent (requestor
supports SCU)
(5528) 06-26 13:20:46.32 MC3 T3: | (73)SCP Role: Not sent (requestor
supports SCP)
(5528) 06-26 13:20:46.32 MC3 T3: | (73)No extended application info
(5528) 06-26 13:20:46.32 MC3 T3: | (74)Pres_context_id: 147
(5528) 06-26 13:20:46.32 MC3 T3: | (74)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.2 (STUDY_ROOT_QR_MOVE)
(5528) 06-26 13:20:46.32 MC3 T3: | (74)Proposed Transfer Syntaxes:
```

```
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.4.51 (JPEG Extended 2
      4)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.4.70 (JPEG Lossless
      Hier 14)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2 (Implicit Little
      Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.1 (Explicit Little
      Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | 1.2.840.10008.1.2.2 (Explicit Big
      Endian)
(5528) 06-26 13:20:46.32 MC3 T3: | (74)SCU Role: Not sent (requestor
      supports SCU)
(5528) 06-26 13:20:46.32 MC3 T3: | (74)SCP Role: Not sent (requestor
      supports SCP)
(5528) 06-26 13:20:46.32 MC3 T3: | (74)No extended application info
```

THE RECEIVING SYSTEM RESPONDS.

```
(5528) 06-26 13:20:46.34 MC3 T3: RESPONSE Association Parameters:
(5528) 06-26 13:20:46.34 MC3 T3: | Calling_presentation_address: ESOFTB
(5528) 06-26 13:20:46.34 MC3 T3: | Called_presentation_address: aqnet_ae:105
(5528) 06-26 13:20:46.34 MC3 T3: | Called_application_title: AQNET_AE
(5528) 06-26 13:20:46.34 MC3 T3: | Calling_application_title: ESOFTB
```

HERE IS THE STANDARD ACN

```
(5528) 06-26 13:20:46.34 MC3 T3: | Application_context_name:
      1.2.840.10008.3.1.1.1
```

MAX PROTOCOL DATA UNIT FOR THE RECEIVER IS 129940

```
(5528) 06-26 13:20:46.34 MC3 T3: | PDU_maximum_length: 129940
```

THIS IS PACS IMPLEMENTATION

```
(5528) 06-26 13:20:46.34 MC3 T3: | Implem_class_uid:
      2.16.840.1.113669.632.21.1
(5528) 06-26 13:20:46.34 MC3 T3: | Version_name: AQNET430-370IB5
(5528) 06-26 13:20:46.34 MC3 T3: | Max_operations_invoked: Not Sent
(5528) 06-26 13:20:46.34 MC3 T3: | Max_operations_performed: Not Sent
```

NOW THE RECEIVER TELLS THE SENDER WHICH OF THE OFFERED

"TRANSFER SYNTAXES" THE RECEIVER WILL ACCEPT. NM IS CONTEXT #13.

```
(5528) 06-26 13:20:46.34 MC3 T3: | PRESENTATION CONTEXTS:
(5528) 06-26 13:20:46.34 MC3 T3: | (1)Pres_context_id: 1
(5528) 06-26 13:20:46.34 MC3 T3: | (1)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1 (STANDARD_CR)
(5528) 06-26 13:20:46.34 MC3 T3: | (1)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.34 MC3 T3: | (1)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (1)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (1)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (1)Negotiation Result: Accepted
(5528) 06-26 13:20:46.34 MC3 T3: | (2)Pres_context_id: 3
(5528) 06-26 13:20:46.34 MC3 T3: | (2)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1 (STANDARD_CR)
(5528) 06-26 13:20:46.34 MC3 T3: | (2)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.34 MC3 T3: | (2)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (2)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (2)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (2)Negotiation Result: Accepted
(5528) 06-26 13:20:46.34 MC3 T3: | (3)Pres_context_id: 5
(5528) 06-26 13:20:46.34 MC3 T3: | (3)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.2 (STANDARD_CT)
(5528) 06-26 13:20:46.34 MC3 T3: | (3)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.34 MC3 T3: | (3)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (3)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (3)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (3)Negotiation Result: Accepted
(5528) 06-26 13:20:46.34 MC3 T3: | (4)Pres_context_id: 7
(5528) 06-26 13:20:46.34 MC3 T3: | (4)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.2 (STANDARD_CT)
(5528) 06-26 13:20:46.34 MC3 T3: | (4)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.34 MC3 T3: | (4)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (4)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (4)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (4)Negotiation Result: Accepted
(5528) 06-26 13:20:46.34 MC3 T3: | (5)Pres_context_id: 9
(5528) 06-26 13:20:46.34 MC3 T3: | (5)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.4 (STANDARD_MR)
(5528) 06-26 13:20:46.34 MC3 T3: | (5)Result_transfer_syntax:
```

```
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.34 MC3 T3: | (5)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (5)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (5)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (5)Negotiation Result: Accepted
(5528) 06-26 13:20:46.34 MC3 T3: | (6)Pres_context_id: 11
(5528) 06-26 13:20:46.34 MC3 T3: | (6)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.4 (STANDARD_MR)
(5528) 06-26 13:20:46.34 MC3 T3: | (6)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:46.34 MC3 T3: | (6)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (6)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (6)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (6)Negotiation Result: Accepted
```

FOR NM, THE RECEIVER ACCEPTS IMPLICIT LITTLE ENDIAN

```
(5528) 06-26 13:20:46.34 MC3 T3: | (7)Pres_context_id: 13
(5528) 06-26 13:20:46.34 MC3 T3: | (7)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.20 (STANDARD_NM)
(5528) 06-26 13:20:46.34 MC3 T3: | (7)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.34 MC3 T3: | (7)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:46.34 MC3 T3: | (7)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:46.34 MC3 T3: | (7)No extended application info
(5528) 06-26 13:20:46.34 MC3 T3: | (7)Negotiation Result: Accepted
(5528) 06-26 13:20:46.34 MC3 T3: | (8)Pres_context_id: 15
(5528) 06-26 13:20:46.34 MC3 T3: | (8)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.20 (STANDARD_NM)
(5528) 06-26 13:20:46.34 MC3 T3: | (8)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:46.34 MC3 T3: | (8)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:46.35 MC3 T3: | (8)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.15 MC3 T3: | (8)No extended application info
(5528) 06-26 13:20:47.15 MC3 T3: | (8)Negotiation Result: Accepted
(5528) 06-26 13:20:47.15 MC3 T3: | (9)Pres_context_id: 17
(5528) 06-26 13:20:47.15 MC3 T3: | (9)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.7 (STANDARD_SEC_CAPTURE)
(5528) 06-26 13:20:47.15 MC3 T3: | (9)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.15 MC3 T3: | (9)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.15 MC3 T3: | (9)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.15 MC3 T3: | (9)No extended application info
(5528) 06-26 13:20:47.15 MC3 T3: | (9)Negotiation Result: Accepted
(5528) 06-26 13:20:47.15 MC3 T3: | (10)Pres_context_id: 19
```

(5528) 06-26 13:20:47.15 MC3 T3: | (10)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.7 (STANDARD_SEC_CAPTURE)
(5528) 06-26 13:20:47.15 MC3 T3: | (10)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.17 MC3 T3: | (10)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (10)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (10)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (10)Negotiation Result: Accepted
(5528) 06-26 13:20:47.17 MC3 T3: | (11)Pres_context_id: 21
(5528) 06-26 13:20:47.17 MC3 T3: | (11)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.6.1 (STANDARD_US)
(5528) 06-26 13:20:47.17 MC3 T3: | (11)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.17 MC3 T3: | (11)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (11)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (11)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (11)Negotiation Result: Accepted
(5528) 06-26 13:20:47.17 MC3 T3: | (12)Pres_context_id: 23
(5528) 06-26 13:20:47.17 MC3 T3: | (12)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.6.1 (STANDARD_US)
(5528) 06-26 13:20:47.17 MC3 T3: | (12)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.17 MC3 T3: | (12)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (12)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (12)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (12)Negotiation Result: Accepted
(5528) 06-26 13:20:47.17 MC3 T3: | (13)Pres_context_id: 25
(5528) 06-26 13:20:47.17 MC3 T3: | (13)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.3.1 (STANDARD_US_MF)
(5528) 06-26 13:20:47.17 MC3 T3: | (13)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.17 MC3 T3: | (13)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (13)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (13)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (13)Negotiation Result: Accepted
(5528) 06-26 13:20:47.17 MC3 T3: | (14)Pres_context_id: 27
(5528) 06-26 13:20:47.17 MC3 T3: | (14)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.3.1 (STANDARD_US_MF)
(5528) 06-26 13:20:47.17 MC3 T3: | (14)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.17 MC3 T3: | (14)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (14)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (14)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (14)Negotiation Result: Accepted
(5528) 06-26 13:20:47.17 MC3 T3: | (15)Pres_context_id: 29
(5528) 06-26 13:20:47.17 MC3 T3: | (15)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.128 (STANDARD_PET)

(5528) 06-26 13:20:47.17 MC3 T3: | (15)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.17 MC3 T3: | (15)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (15)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (15)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (15)Negotiation Result: Accepted
(5528) 06-26 13:20:47.17 MC3 T3: | (16)Pres_context_id: 31
(5528) 06-26 13:20:47.17 MC3 T3: | (16)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.128 (STANDARD_PET)
(5528) 06-26 13:20:47.17 MC3 T3: | (16)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.17 MC3 T3: | (16)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.17 MC3 T3: | (16)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.17 MC3 T3: | (16)No extended application info
(5528) 06-26 13:20:47.17 MC3 T3: | (16)Negotiation Result: Accepted

SOME STUDY TYPES ARE REJECTED

(5528) 06-26 13:20:47.17 MC3 T3: | (17)Pres_context_id: 33
(5528) 06-26 13:20:47.17 MC3 T3: | (17)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.1 (STANDARD_WAVEFORM_12_LEAD_ECG)
(5528) 06-26 13:20:47.17 MC3 T3: | (17)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.17 MC3 T3: | (18)Pres_context_id: 35
(5528) 06-26 13:20:47.17 MC3 T3: | (18)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.1 (STANDARD_WAVEFORM_12_LEAD_ECG)
(5528) 06-26 13:20:47.17 MC3 T3: | (18)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.17 MC3 T3: | (19)Pres_context_id: 37
(5528) 06-26 13:20:47.17 MC3 T3: | (19)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.2 (STANDARD_WAVEFORM_GENERAL_ECG)
(5528) 06-26 13:20:47.17 MC3 T3: | (19)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.17 MC3 T3: | (20)Pres_context_id: 39
(5528) 06-26 13:20:47.17 MC3 T3: | (20)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.2 (STANDARD_WAVEFORM_GENERAL_ECG)
(5528) 06-26 13:20:47.17 MC3 T3: | (20)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.17 MC3 T3: | (21)Pres_context_id: 41
(5528) 06-26 13:20:47.17 MC3 T3: | (21)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.3
(STANDARD_WAVEFORM_AMBULATORY_ECG)
(5528) 06-26 13:20:47.17 MC3 T3: | (21)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.17 MC3 T3: | (22)Pres_context_id: 43
(5528) 06-26 13:20:47.17 MC3 T3: | (22)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.1.3
(STANDARD_WAVEFORM_AMBULATORY_ECG)
(5528) 06-26 13:20:47.18 MC3 T3: | (22)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (23)Pres_context_id: 45
(5528) 06-26 13:20:47.18 MC3 T3: | (23)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.9.2.1 (STANDARD_WAVEFORM_HEMODYNAMIC)


```
(5528) 06-26 13:20:47.18 MC3 T3: | (23)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (24)Pres_context_id: 47
(5528) 06-26 13:20:47.18 MC3 T3: | (24)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.2.1 (STANDARD_WAVEFORM_HEMODYNAMIC)
(5528) 06-26 13:20:47.18 MC3 T3: | (24)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (25)Pres_context_id: 49
(5528) 06-26 13:20:47.18 MC3 T3: | (25)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.3.1 (STANDARD_WAVEFORM_CARDIAC_EP)
(5528) 06-26 13:20:47.18 MC3 T3: | (25)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (26)Pres_context_id: 51
(5528) 06-26 13:20:47.18 MC3 T3: | (26)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.3.1 (STANDARD_WAVEFORM_CARDIAC_EP)
(5528) 06-26 13:20:47.18 MC3 T3: | (26)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (27)Pres_context_id: 53
(5528) 06-26 13:20:47.18 MC3 T3: | (27)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.4.1
      (STANDARD_WAVEFORM_BASIC_VOICE_AU)
(5528) 06-26 13:20:47.18 MC3 T3: | (27)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (28)Pres_context_id: 55
(5528) 06-26 13:20:47.18 MC3 T3: | (28)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.9.4.1
      (STANDARD_WAVEFORM_BASIC_VOICE_AU)
(5528) 06-26 13:20:47.18 MC3 T3: | (28)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.18 MC3 T3: | (29)Pres_context_id: 57
(5528) 06-26 13:20:47.18 MC3 T3: | (29)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.12.1 (STANDARD_XRAY_ANGIO)
(5528) 06-26 13:20:47.18 MC3 T3: | (29)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.18 MC3 T3: | (29)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.18 MC3 T3: | (29)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.18 MC3 T3: | (29)No extended application info
(5528) 06-26 13:20:47.18 MC3 T3: | (29)Negotiation Result: Accepted
(5528) 06-26 13:20:47.18 MC3 T3: | (30)Pres_context_id: 59
(5528) 06-26 13:20:47.18 MC3 T3: | (30)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.12.1 (STANDARD_XRAY_ANGIO)
(5528) 06-26 13:20:47.18 MC3 T3: | (30)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.18 MC3 T3: | (30)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.18 MC3 T3: | (30)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.18 MC3 T3: | (30)No extended application info
(5528) 06-26 13:20:47.18 MC3 T3: | (30)Negotiation Result: Accepted
(5528) 06-26 13:20:47.18 MC3 T3: | (31)Pres_context_id: 61
(5528) 06-26 13:20:47.18 MC3 T3: | (31)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.12.3 (STANDARD_XRAY_ANGIO_BIPLANE)
(5528) 06-26 13:20:47.18 MC3 T3: | (31)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.18 MC3 T3: | (31)SCU Role: Not sent (remote accepts
```

```
requestor SCU role)
(5528) 06-26 13:20:47.18 MC3 T3: | (31)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.18 MC3 T3: | (31)No extended application info
(5528) 06-26 13:20:47.18 MC3 T3: | (31)Negotiation Result: Accepted
(5528) 06-26 13:20:47.18 MC3 T3: | (32)Pres_context_id: 63
(5528) 06-26 13:20:47.18 MC3 T3: | (32)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.3 (STANDARD_XRAY_ANGIO_BIPLANE)
(5528) 06-26 13:20:47.18 MC3 T3: | (32)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.18 MC3 T3: | (32)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.18 MC3 T3: | (32)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.18 MC3 T3: | (32)No extended application info
(5528) 06-26 13:20:47.18 MC3 T3: | (32)Negotiation Result: Accepted
(5528) 06-26 13:20:47.18 MC3 T3: | (33)Pres_context_id: 65
(5528) 06-26 13:20:47.18 MC3 T3: | (33)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.2 (STANDARD_XRAY_RF)
(5528) 06-26 13:20:47.18 MC3 T3: | (33)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.18 MC3 T3: | (33)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.18 MC3 T3: | (33)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.18 MC3 T3: | (33)No extended application info
(5528) 06-26 13:20:47.18 MC3 T3: | (33)Negotiation Result: Accepted
(5528) 06-26 13:20:47.18 MC3 T3: | (34)Pres_context_id: 67
(5528) 06-26 13:20:47.18 MC3 T3: | (34)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.12.2 (STANDARD_XRAY_RF)
(5528) 06-26 13:20:47.18 MC3 T3: | (34)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.18 MC3 T3: | (34)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.18 MC3 T3: | (34)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.18 MC3 T3: | (34)No extended application info
(5528) 06-26 13:20:47.18 MC3 T3: | (34)Negotiation Result: Accepted
(5528) 06-26 13:20:47.20 MC3 T3: | (35)Pres_context_id: 69
(5528) 06-26 13:20:47.20 MC3 T3: | (35)Abstract_syntax: 1.3.12.2.1107.5.9.1
(PRIVATE_MEDCOM_CSA_NON_IMAGE)
(5528) 06-26 13:20:47.20 MC3 T3: | (35)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (36)Pres_context_id: 71
(5528) 06-26 13:20:47.20 MC3 T3: | (36)Abstract_syntax: 1.3.12.2.1107.5.9.1
(PRIVATE_MEDCOM_CSA_NON_IMAGE)
(5528) 06-26 13:20:47.20 MC3 T3: | (36)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (37)Pres_context_id: 73
(5528) 06-26 13:20:47.20 MC3 T3: | (37)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.3 (STANDARD_RT_STRUCTURE_SET)
(5528) 06-26 13:20:47.20 MC3 T3: | (37)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (38)Pres_context_id: 75
(5528) 06-26 13:20:47.20 MC3 T3: | (38)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.3 (STANDARD_RT_STRUCTURE_SET)
(5528) 06-26 13:20:47.20 MC3 T3: | (38)Negotiation Result: REJECTED:
```

```
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (39)Pres_context_id: 77
(5528) 06-26 13:20:47.20 MC3 T3: | (39)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1.1 (STANDARD_DX_PRESENT)
(5528) 06-26 13:20:47.20 MC3 T3: | (39)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.20 MC3 T3: | (39)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.20 MC3 T3: | (39)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.20 MC3 T3: | (39)No extended application info
(5528) 06-26 13:20:47.20 MC3 T3: | (39)Negotiation Result: Accepted
(5528) 06-26 13:20:47.20 MC3 T3: | (40)Pres_context_id: 79
(5528) 06-26 13:20:47.20 MC3 T3: | (40)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1.1 (STANDARD_DX_PRESENT)
(5528) 06-26 13:20:47.20 MC3 T3: | (40)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.20 MC3 T3: | (40)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.20 MC3 T3: | (40)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.20 MC3 T3: | (40)No extended application info
(5528) 06-26 13:20:47.20 MC3 T3: | (40)Negotiation Result: Accepted
(5528) 06-26 13:20:47.20 MC3 T3: | (41)Pres_context_id: 81
(5528) 06-26 13:20:47.20 MC3 T3: | (41)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1.1.1 (STANDARD_DX_PROCESS)
(5528) 06-26 13:20:47.20 MC3 T3: | (41)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.20 MC3 T3: | (41)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.20 MC3 T3: | (41)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.20 MC3 T3: | (41)No extended application info
(5528) 06-26 13:20:47.20 MC3 T3: | (41)Negotiation Result: Accepted
(5528) 06-26 13:20:47.20 MC3 T3: | (42)Pres_context_id: 83
(5528) 06-26 13:20:47.20 MC3 T3: | (42)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1.1.1 (STANDARD_DX_PROCESS)
(5528) 06-26 13:20:47.20 MC3 T3: | (42)Result_transfer_syntax:
      1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.20 MC3 T3: | (42)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.20 MC3 T3: | (42)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.20 MC3 T3: | (42)No extended application info
(5528) 06-26 13:20:47.20 MC3 T3: | (42)Negotiation Result: Accepted
(5528) 06-26 13:20:47.20 MC3 T3: | (43)Pres_context_id: 85
(5528) 06-26 13:20:47.20 MC3 T3: | (43)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1.1.3 (STANDARD_IO_PRESENT)
(5528) 06-26 13:20:47.20 MC3 T3: | (43)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (44)Pres_context_id: 87
(5528) 06-26 13:20:47.20 MC3 T3: | (44)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.1.1.3 (STANDARD_IO_PRESENT)
(5528) 06-26 13:20:47.20 MC3 T3: | (44)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (45)Pres_context_id: 89
(5528) 06-26 13:20:47.20 MC3 T3: | (45)Abstract_syntax:
```

```
1.2.840.10008.5.1.4.1.1.1.3.1 (STANDARD_IO_PROCESS)
(5528) 06-26 13:20:47.20 MC3 T3: | (45)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (46)Pres_context_id: 91
(5528) 06-26 13:20:47.20 MC3 T3: | (46)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.3.1 (STANDARD_IO_PROCESS)
(5528) 06-26 13:20:47.20 MC3 T3: | (46)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.20 MC3 T3: | (47)Pres_context_id: 93
(5528) 06-26 13:20:47.20 MC3 T3: | (47)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2 (STANDARD_MG_PRESENT)
(5528) 06-26 13:20:47.20 MC3 T3: | (47)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.20 MC3 T3: | (47)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.20 MC3 T3: | (47)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.20 MC3 T3: | (47)No extended application info
(5528) 06-26 13:20:47.20 MC3 T3: | (47)Negotiation Result: Accepted
(5528) 06-26 13:20:47.20 MC3 T3: | (48)Pres_context_id: 95
(5528) 06-26 13:20:47.20 MC3 T3: | (48)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2 (STANDARD_MG_PRESENT)
(5528) 06-26 13:20:47.20 MC3 T3: | (48)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.21 MC3 T3: | (48)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.21 MC3 T3: | (48)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.21 MC3 T3: | (48)No extended application info
(5528) 06-26 13:20:47.21 MC3 T3: | (48)Negotiation Result: Accepted
(5528) 06-26 13:20:47.21 MC3 T3: | (49)Pres_context_id: 97
(5528) 06-26 13:20:47.21 MC3 T3: | (49)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2.1 (STANDARD_MG_PROCESS)
(5528) 06-26 13:20:47.21 MC3 T3: | (49)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.21 MC3 T3: | (49)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.21 MC3 T3: | (49)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.21 MC3 T3: | (49)No extended application info
(5528) 06-26 13:20:47.21 MC3 T3: | (49)Negotiation Result: Accepted
(5528) 06-26 13:20:47.21 MC3 T3: | (50)Pres_context_id: 99
(5528) 06-26 13:20:47.21 MC3 T3: | (50)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.1.2.1 (STANDARD_MG_PROCESS)
(5528) 06-26 13:20:47.21 MC3 T3: | (50)Result_transfer_syntax:
1.2.840.10008.1.2.4.70 (JPEG Lossless Hier 14)
(5528) 06-26 13:20:47.21 MC3 T3: | (50)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.21 MC3 T3: | (50)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.21 MC3 T3: | (50)No extended application info
(5528) 06-26 13:20:47.21 MC3 T3: | (50)Negotiation Result: Accepted
(5528) 06-26 13:20:47.21 MC3 T3: | (51)Pres_context_id: 101
(5528) 06-26 13:20:47.21 MC3 T3: | (51)Abstract_syntax:
1.2.840.10008.5.1.4.1.1.481.2 (STANDARD_RT_DOSE)
(5528) 06-26 13:20:47.21 MC3 T3: | (51)Negotiation Result: REJECTED:
Abstract syntax
```

```
(5528) 06-26 13:20:47.21 MC3 T3: | (52)Pres_context_id: 103
(5528) 06-26 13:20:47.21 MC3 T3: | (52)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.481.2 (STANDARD_RT_DOSE)
(5528) 06-26 13:20:47.21 MC3 T3: | (52)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.21 MC3 T3: | (53)Pres_context_id: 105
(5528) 06-26 13:20:47.21 MC3 T3: | (53)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.481.1 (STANDARD_RT_IMAGE)
(5528) 06-26 13:20:47.21 MC3 T3: | (53)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.21 MC3 T3: | (54)Pres_context_id: 107
(5528) 06-26 13:20:47.21 MC3 T3: | (54)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.481.1 (STANDARD_RT_IMAGE)
(5528) 06-26 13:20:47.21 MC3 T3: | (54)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.21 MC3 T3: | (55)Pres_context_id: 109
(5528) 06-26 13:20:47.21 MC3 T3: | (55)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.481.5 (STANDARD_RT_PLAN)
(5528) 06-26 13:20:47.21 MC3 T3: | (55)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.21 MC3 T3: | (56)Pres_context_id: 111
(5528) 06-26 13:20:47.21 MC3 T3: | (56)Abstract_syntax:
      1.2.840.10008.5.1.4.1.1.481.5 (STANDARD_RT_PLAN)
(5528) 06-26 13:20:47.21 MC3 T3: | (56)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.21 MC3 T3: | (57)Pres_context_id: 113
(5528) 06-26 13:20:47.21 MC3 T3: | (57)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.1 (PATIENT_ROOT_QR_FIND)
(5528) 06-26 13:20:47.21 MC3 T3: | (57)Result_transfer_syntax:
      1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.21 MC3 T3: | (57)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.21 MC3 T3: | (57)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.21 MC3 T3: | (57)No extended application info
(5528) 06-26 13:20:47.21 MC3 T3: | (57)Negotiation Result: Accepted
(5528) 06-26 13:20:47.21 MC3 T3: | (58)Pres_context_id: 115
(5528) 06-26 13:20:47.21 MC3 T3: | (58)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.1 (PATIENT_ROOT_QR_FIND)
(5528) 06-26 13:20:47.21 MC3 T3: | (58)Result_transfer_syntax:
      1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.21 MC3 T3: | (58)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.21 MC3 T3: | (58)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.21 MC3 T3: | (58)No extended application info
(5528) 06-26 13:20:47.21 MC3 T3: | (58)Negotiation Result: Accepted
(5528) 06-26 13:20:47.21 MC3 T3: | (59)Pres_context_id: 117
(5528) 06-26 13:20:47.21 MC3 T3: | (59)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.2.1 (STUDY_ROOT_QR_FIND)
(5528) 06-26 13:20:47.21 MC3 T3: | (59)Result_transfer_syntax:
      1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.21 MC3 T3: | (59)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.21 MC3 T3: | (59)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.21 MC3 T3: | (59)No extended application info
```

```
(5528) 06-26 13:20:47.21 MC3 T3: | (59)Negotiation Result: Accepted
(5528) 06-26 13:20:47.21 MC3 T3: | (60)Pres_context_id: 119
(5528) 06-26 13:20:47.21 MC3 T3: | (60)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.2.1 (STUDY_ROOT_QR_FIND)
(5528) 06-26 13:20:47.21 MC3 T3: | (60)Result_transfer_syntax:
      1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.21 MC3 T3: | (60)SCU Role: Not sent (remote accepts
      requestor SCU role)
(5528) 06-26 13:20:47.23 MC3 T3: | (60)SCP Role: Not sent (remote accepts
      requestor SCP role)
(5528) 06-26 13:20:47.23 MC3 T3: | (60)No extended application info
(5528) 06-26 13:20:47.23 MC3 T3: | (60)Negotiation Result: Accepted
(5528) 06-26 13:20:47.23 MC3 T3: | (61)Pres_context_id: 121
(5528) 06-26 13:20:47.23 MC3 T3: | (61)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.3.1 (PATIENT_STUDY_ONLY_QR_FIND)
(5528) 06-26 13:20:47.23 MC3 T3: | (61)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (62)Pres_context_id: 123
(5528) 06-26 13:20:47.23 MC3 T3: | (62)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.3.1 (PATIENT_STUDY_ONLY_QR_FIND)
(5528) 06-26 13:20:47.23 MC3 T3: | (62)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (63)Pres_context_id: 125
(5528) 06-26 13:20:47.23 MC3 T3: | (63)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.3 (PATIENT_ROOT_QR_GET)
(5528) 06-26 13:20:47.23 MC3 T3: | (63)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (64)Pres_context_id: 127
(5528) 06-26 13:20:47.23 MC3 T3: | (64)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.3 (PATIENT_ROOT_QR_GET)
(5528) 06-26 13:20:47.23 MC3 T3: | (64)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (65)Pres_context_id: 129
(5528) 06-26 13:20:47.23 MC3 T3: | (65)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.3.3 (PATIENT_STUDY_ONLY_QR_GET)
(5528) 06-26 13:20:47.23 MC3 T3: | (65)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (66)Pres_context_id: 131
(5528) 06-26 13:20:47.23 MC3 T3: | (66)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.3.3 (PATIENT_STUDY_ONLY_QR_GET)
(5528) 06-26 13:20:47.23 MC3 T3: | (66)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (67)Pres_context_id: 133
(5528) 06-26 13:20:47.23 MC3 T3: | (67)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.2.3 (STUDY_ROOT_QR_GET)
(5528) 06-26 13:20:47.23 MC3 T3: | (67)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (68)Pres_context_id: 135
(5528) 06-26 13:20:47.23 MC3 T3: | (68)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.2.3 (STUDY_ROOT_QR_GET)
(5528) 06-26 13:20:47.23 MC3 T3: | (68)Negotiation Result: REJECTED:
      Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (69)Pres_context_id: 137
(5528) 06-26 13:20:47.23 MC3 T3: | (69)Abstract_syntax:
      1.2.840.10008.5.1.4.1.2.1.2 (PATIENT_ROOT_QR_MOVE)
(5528) 06-26 13:20:47.23 MC3 T3: | (69)Result_transfer_syntax:
      1.2.840.10008.1.2 (Implicit Little Endian)
```

```
(5528) 06-26 13:20:47.23 MC3 T3: | (69)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.23 MC3 T3: | (69)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.23 MC3 T3: | (69)No extended application info
(5528) 06-26 13:20:47.23 MC3 T3: | (69)Negotiation Result: Accepted
(5528) 06-26 13:20:47.23 MC3 T3: | (70)Pres_context_id: 139
(5528) 06-26 13:20:47.23 MC3 T3: | (70)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.1.2 (PATIENT_ROOT_QR_MOVE)
(5528) 06-26 13:20:47.23 MC3 T3: | (70)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.23 MC3 T3: | (70)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.23 MC3 T3: | (70)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.23 MC3 T3: | (70)No extended application info
(5528) 06-26 13:20:47.23 MC3 T3: | (70)Negotiation Result: Accepted
(5528) 06-26 13:20:47.23 MC3 T3: | (71)Pres_context_id: 141
(5528) 06-26 13:20:47.23 MC3 T3: | (71)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.3.2 (PATIENT_STUDY_ONLY_QR_MOVE)
(5528) 06-26 13:20:47.23 MC3 T3: | (71)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (72)Pres_context_id: 143
(5528) 06-26 13:20:47.23 MC3 T3: | (72)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.3.2 (PATIENT_STUDY_ONLY_QR_MOVE)
(5528) 06-26 13:20:47.23 MC3 T3: | (72)Negotiation Result: REJECTED:
Abstract syntax
(5528) 06-26 13:20:47.23 MC3 T3: | (73)Pres_context_id: 145
(5528) 06-26 13:20:47.23 MC3 T3: | (73)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.2 (STUDY_ROOT_QR_MOVE)
(5528) 06-26 13:20:47.23 MC3 T3: | (73)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.23 MC3 T3: | (73)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.23 MC3 T3: | (73)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.23 MC3 T3: | (73)No extended application info
(5528) 06-26 13:20:47.23 MC3 T3: | (73)Negotiation Result: Accepted
(5528) 06-26 13:20:47.23 MC3 T3: | (74)Pres_context_id: 147
(5528) 06-26 13:20:47.23 MC3 T3: | (74)Abstract_syntax:
1.2.840.10008.5.1.4.1.2.2.2 (STUDY_ROOT_QR_MOVE)
(5528) 06-26 13:20:47.23 MC3 T3: | (74)Result_transfer_syntax:
1.2.840.10008.1.2 (Implicit Little Endian)
(5528) 06-26 13:20:47.23 MC3 T3: | (74)SCU Role: Not sent (remote accepts
requestor SCU role)
(5528) 06-26 13:20:47.23 MC3 T3: | (74)SCP Role: Not sent (remote accepts
requestor SCP role)
(5528) 06-26 13:20:47.23 MC3 T3: | (74)No extended application info
(5528) 06-26 13:20:47.23 MC3 T3: | (74)Negotiation Result: Accepted
```

REJECTIONS:

ENSURE NM WAS NOT REJECTED

NM IS CONTEXT # 13 & #15

(5528) 06-26 13:20:47.25 MC3 W: Presentation context 33 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.1.1 (STANDARD_WAVEFORM_12_LEAD_ECG)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 35 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.1.1 (STANDARD_WAVEFORM_12_LEAD_ECG)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 37 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.1.2 (STANDARD_WAVEFORM_GENERAL_ECG)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 39 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.1.2 (STANDARD_WAVEFORM_GENERAL_ECG)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 41 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.1.3
(STANDARD_WAVEFORM_AMBULATORY_ECG)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 43 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.1.3
(STANDARD_WAVEFORM_AMBULATORY_ECG)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 45 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.2.1 (STANDARD_WAVEFORM_HEMODYNAMIC)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 47 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.2.1 (STANDARD_WAVEFORM_HEMODYNAMIC)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 49 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.3.1 (STANDARD_WAVEFORM_CARDIAC_EP)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 51 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.3.1 (STANDARD_WAVEFORM_CARDIAC_EP)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 53 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.4.1
(STANDARD_WAVEFORM_BASIC_VOICE_AU)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 55 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.9.4.1
(STANDARD_WAVEFORM_BASIC_VOICE_AU)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 69 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax: 1.3.12.2.1107.5.9.1
(PRIVATE_MEDCOM_CSA_NON_IMAGE)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 71 rejected, reason 0x03

(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax: 1.3.12.2.1107.5.9.1
(PRIVATE_MEDCOM_CSA_NON_IMAGE)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 73 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.3 (STANDARD_RT_STRUCTURE_SET)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 75 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.3 (STANDARD_RT_STRUCTURE_SET)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 85 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.1.3 (STANDARD_IO_PRESENT)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 87 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.1.3 (STANDARD_IO_PRESENT)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 89 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.1.3.1 (STANDARD_IO_PROCESS)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 91 rejected, reason 0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.1.3.1 (STANDARD_IO_PROCESS)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 101 rejected, reason
0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.2 (STANDARD_RT_DOSE)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 103 rejected, reason
0x03
(5528) 06-26 13:20:47.25 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.2 (STANDARD_RT_DOSE)
(5528) 06-26 13:20:47.25 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.25 MC3 W: Presentation context 105 rejected, reason
0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.1 (STANDARD_RT_IMAGE)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 107 rejected, reason
0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.1 (STANDARD_RT_IMAGE)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 109 rejected, reason
0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.5 (STANDARD_RT_PLAN)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 111 rejected, reason
0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.1.481.5 (STANDARD_RT_PLAN)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected

(5528) 06-26 13:20:47.26 MC3 W: Presentation context 121 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.3.1 (PATIENT_STUDY_ONLY_QR_FIND)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 123 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.3.1 (PATIENT_STUDY_ONLY_QR_FIND)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 125 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.1.3 (PATIENT_ROOT_QR_GET)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 127 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.1.3 (PATIENT_ROOT_QR_GET)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 129 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.3.3 (PATIENT_STUDY_ONLY_QR_GET)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 131 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.3.3 (PATIENT_STUDY_ONLY_QR_GET)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 133 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.2.3 (STUDY_ROOT_QR_GET)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 135 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.2.3 (STUDY_ROOT_QR_GET)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 141 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.3.2 (PATIENT_STUDY_ONLY_QR_MOVE)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected
(5528) 06-26 13:20:47.26 MC3 W: Presentation context 143 rejected, reason 0x03
(5528) 06-26 13:20:47.26 MC3 W: | Abstract syntax:
1.2.840.10008.5.1.4.1.2.3.2 (PATIENT_STUDY_ONLY_QR_MOVE)
(5528) 06-26 13:20:47.26 MC3 W: | Abstract Syntax Rejected

PERFORM A STANDARD NM STORE USING CONTEXT #13

(5528) 06-26 13:20:47.26 MC3 T7: Message (STANDARD_NM,C_STORE_RQ) sent on context 13

(5528) 06-26 13:20:47.26 MC3 T7: Group 0x0000 Elements:
(5528) 06-26 13:20:47.26 MC3 T7: (0000,0000) Group 0000 Length
(5528) 06-26 13:20:47.26 MC3 T7: (1): 136

THIS IS A NM STUDY

(5528) 06-26 13:20:47.26 MC3 T7: (0000,0002) Affected SOP Class UID
(5528) 06-26 13:20:47.26 MC3 T7: (1): 1.2.840.10008.5.1.4.1.1.20
(5528) 06-26 13:20:47.26 MC3 T7: (0000,0100) Command Field
(5528) 06-26 13:20:47.26 MC3 T7: (1): 1
(5528) 06-26 13:20:47.26 MC3 T7: (0000,0110) Message ID
(5528) 06-26 13:20:47.26 MC3 T7: (1): 0
(5528) 06-26 13:20:47.26 MC3 T7: (0000,0700) Priority
(5528) 06-26 13:20:47.26 MC3 T7: (1): 0
(5528) 06-26 13:20:47.26 MC3 T7: (0000,0800) Data Set Type
(5528) 06-26 13:20:47.26 MC3 T7: (1): 0

UNIQUE IMAGE ID

(5528) 06-26 13:20:47.26 MC3 T7: (0000,1000) Affected SOP Instance UID
(5528) 06-26 13:20:47.26 MC3 T7: (1):
1.3.12.2.1107.5.6.1.6535.32510109062617003339000000002
(5528) 06-26 13:20:48.42 MC3 T2: Message (STANDARD_NM,C_STORE_RQ) sent on
context 13

LIST THE DICOM ELEMENTS OF THE TRANSFERRED FILE

```
(5528)
=====
=====
(5528) Message ID: 100
(5528) Service Supported: STANDARD_NM (0061)
(5528) Command Supported: C_STORE_RQ
(5528) Message Attributes:
(5528) 0000,0000 Group 0000 Length VR: UL
      VM: 1
(5528) (0000004) |136|
(5528) 0000,0002 Affected SOP Class UID VR: UI
      VM: 1
(5528) (0000026) |1.2.840.10008.5.1.4.1.1.20|
(5528) 0000,0100 Command Field VR: US
      VM: 1
(5528) (0000002) |1|
(5528) 0000,0110 Message ID VR: US
      VM: 1
(5528) (0000002) |0|
(5528) 0000,0700 Priority VR: US
      VM: 1
(5528) (0000002) |0|
(5528) 0000,0800 Data Set Type VR: US
      VM: 1
(5528) (0000002) |0|
(5528) 0000,1000 Affected SOP Instance UID VR: UI
      VM: 1
(5528) (0000054)
```

```
|1.3.12.2.1107.5.6.1.6535.325101090626170033390000|
(5528) |0002|
(5528) 0008,0005 Specific Character Set VR: CS
      VM: 1-n
(5528) (0000010) |ISO_IR 100|
(5528) 0008,0008 Image Type VR: CS
      VM: 1-n
(5528) (0000027) |ORIGINAL|
(5528) |PRIMARY|
(5528) |TOMO|
(5528) |EMISSION|
(5528) 0008,0016 SOP Class UID VR: UI
      VM: 1
(5528) (0000026) |1.2.840.10008.5.1.4.1.1.20|
(5528) 0008,0018 SOP Instance UID VR: UI
      VM: 1
(5528) (0000054)
      |1.3.12.2.1107.5.6.1.6535.325101090626170033390000|
(5528) |0002|
(5528) 0008,0020 Study Date VR: DA
      VM: 1
(5528) (0000008) |20090626|
(5528) 0008,0021 Series Date VR: DA
      VM: 1
(5528) (0000008) |20090626|
(5528) 0008,0022 Acquisition Date VR: DA
      VM: 1
(5528) (0000008) |20090626|
(5528) 0008,0023 Content Date VR: DA
      VM: 1
(5528) (0000008) |20090626|
(5528) 0008,0030 Study Time VR: TM
      VM: 1
(5528) (0000013) |092123.843000|
(5528) 0008,0031 Series Time VR: TM
      VM: 1
(5528) (0000013) |130346.000000|
(5528) 0008,0032 Acquisition Time VR: TM
      VM: 1
(5528) (0000013) |130346.000000|
(5528) 0008,0033 Content Time VR: TM
      VM: 1
(5528) (0000013) |130346.000000|
(5528) 0008,0050 Accession Number VR: SH
      VM: 1
(5528) (0000000) <null>
(5528) 0008,0060 Modality VR: CS
      VM: 1
(5528) (0000002) |NM|
(5528) 0008,0070 Manufacturer VR: LO
      VM: 1
(5528) (0000010) |SIEMENS NM|
(5528) 0008,0090 Referring Physician's Name VR: PN
      VM: 1
(5528) (0000000) <null>
(5528) 0008,1010 Station Name VR: SH
      VM: 1
```

```
(5528)          (0000006) |ESOFTB|
(5528) 0008,1030 Study Description                      VR: LO
                VM: 1
(5528)          (0000020) |Myocardial Perfusion|
(5528) 0008,103E Series Description                    VR: LO
                VM: 1
(5528)          (0000004) |Rest|
(5528) 0008,1070 Operators' Name                       VR: PN
                VM: 1-n
(5528)          (0000000) <null>
(5528) 0008,1090 Manufacturer's Model Name            VR: LO
                VM: 1
(5528)          (0000003) |IP2|
(5528) 0008,2218 Anatomic Region Sequence             VR: SQ
                VM: 1
(5528)
=====
(5528) <Sequence: Item 104 of Message 100 (0008,2218)>
(5528)
=====
(5528) Item ID: 104
(5528) Item Attributes:
(5528) 0008,0100 Code Value                            VR:
                SH VM: 1
                (0000007) |T-32000|
(5528) 0008,0102 Coding Scheme Designator            VR:
                SH VM: 1
                (0000005) |99SDM|
(5528) 0008,0104 Code Meaning                        VR:
                LO VM: 1
                (0000007) |Cardiac|
(5528) ===== End of Item ID: 104 =====
(5528) Message ID: 100 (continued)
(5528) 0010,0010 Patient's Name                      VR: PN
                VM: 1
                (0000010) |SMITH|JOHN|
(5528) 0010,0020 Patient ID                          VR: LO
                VM: 1
                (0000008) |NM012345|
(5528) 0010,0030 Patient's Birth Date                VR: DA
                VM: 1
                (0000008) |19490912|
(5528) 0010,0040 Patient's Sex                      VR: CS
                VM: 1
                (0000001) |M|
(5528) 0010,1010 Patient's Age                      VR: AS
                VM: 1
                (0000004) |059Y|
(5528) 0018,0026 Intervention Drug Information Sequence VR: SQ
                VM: 1
(5528)
=====
(5528) <Sequence: Item 110 of Message 100 (0018,0026)>
(5528)
```

```
=====
=====
(5528) Item ID: 110
(5528) Item Attributes:
(5528) 0018,0028 Intervention Drug Dose VR:
      DS VM: 1
(5528) (0000001) |0|
(5528) 0018,0034 Intervention Drug Name VR:
      LO VM: 1
(5528) (0000000) <null>
(5528) 0018,0035 Intervention Drug Start Time VR:
      TM VM: 1
(5528) (0000013) |000000.000000|
(5528) ===== End of Item ID: 110 =====
(5528) Message ID: 100 (continued)
(5528) 0018,0070 Counts Accumulated VR: IS
      VM: 1
(5528) (0000000) <null>
(5528) 0018,0071 Acquisition Termination Condition VR: CS
      VM: 1
(5528) (0000004) |TIME|
(5528) 0018,1020 Software Version(s) VR: LO
      VM: 1-n
(5528) (0000014) |3.5.7.13|
(5528) (0000014) |5.5.12|
(5528) 0018,1030 Protocol Name VR: LO
      VM: 1
(5528) (0000014) |25108.0.46211@|
(5528) 0018,1200 Date of Last Calibration VR: DA
      VM: 1-n
(5528) (0000048) |20090623|
(5528) (0000048) |20090622|
(5528) (0000048) |20090622|
(5528) (0000048) |20090623|
(5528) (0000048) |20090622|
(5528) (0000048) |20090622|
(5528) 0018,1201 Time of Last Calibration VR: TM
      VM: 1-n
(5528) (0000078) |171237.000000|
(5528) (0000078) |115522.000000|
(5528) (0000078) |123123.000000|
(5528) (0000078) |171237.000000|
(5528) (0000078) |115522.000000|
(5528) (0000078) |123123.000000|
(5528) 0018,5100 Patient Position VR: CS
      VM: 1
(5528) (0000003) |FFS|
(5528) 0019,0010 Private creator code VR: LO
      VM: 1
(5528) (0000014) |SIEMENS MED NM|
(5528) 0019,100F Siemens ICON Data Type RETIRED VR: SL
      VM: 1-n
(5528) (0000004) |8|
(5528) 0020,000D Study Instance UID VR: UI
      VM: 1
(5528) (0000054)
      |1.3.12.2.1107.5.6.1.8813.3250010906261728091710000|
```

```
(5528) |0002|
(5528) 0020,000E Series Instance UID VR: UI
      VM: 1
(5528) (0000054)
      |1.3.12.2.1107.5.6.1.6535.325101090626170033390000|
(5528) |0000|
(5528) 0020,0010 Study ID VR: SH
      VM: 1
(5528) (0000001) |1|
(5528) 0020,0011 Series Number VR: IS
      VM: 1
(5528) (0000001) |1|
(5528) 0020,0013 Instance Number VR: IS
      VM: 1
(5528) (0000001) |1|
(5528) 0020,0052 Frame of Reference UID VR: UI
      VM: 1
(5528) (0000054)
      |1.3.12.2.1107.5.6.1.6535.325101090626170033390000|
(5528) |0001|
(5528) 0028,0002 Samples per Pixel VR: US
      VM: 1
(5528) (0000002) |1|
(5528) 0028,0004 Photometric Interpretation VR: CS
      VM: 1
(5528) (0000011) |MONOCHROME2|
(5528) 0028,0008 Number of Frames VR: IS
      VM: 1
(5528) (0000002) |32|
(5528) 0028,0009 Frame Increment Pointer VR: AT
      VM: 1-n
(5528) (0000016) |5505040|
(5528) |5505056|
(5528) |5505104|
(5528) |5505168|
(5528) 0028,0010 Rows VR: US
      VM: 1
(5528) (0000002) |64|
(5528) 0028,0011 Columns VR: US
      VM: 1
(5528) (0000002) |64|
(5528) 0028,0030 Pixel Spacing VR: DS
      VM: 2
(5528) (0000030) |6.5913400650024|
(5528) |6.5913400650024|
(5528) 0028,0051 Corrected Image VR: CS
      VM: 1-n
(5528) (0000012) |UNIF|
(5528) |UNIF|
(5528) |UNIF|
(5528) 0028,0100 Bits Allocated VR: US
      VM: 1
(5528) (0000002) |16|
(5528) 0028,0101 Bits Stored VR: US
      VM: 1
(5528) (0000002) |16|
(5528) 0028,0102 High Bit VR: US
```

```
VM: 1
(5528) (0000002) |15|
(5528) 0028,0103 Pixel Representation VR: US
VM: 1
(5528) (0000002) |0|
(5528) 0028,0106 Smallest Image Pixel Value VR: US
VM: 1
(5528) (0000002) |0|
(5528) 0028,0107 Largest Image Pixel Value VR: US
VM: 1
(5528) (0000002) |112|
(5528) 0028,1050 Window Center VR: DS
VM: 1-n
(5528) (0000002) |56|
(5528) 0028,1051 Window Width VR: DS
VM: 1-n
(5528) (0000003) |112|
(5528) 0029,0010 Private creator code VR: LO
VM: 1
(5528) (0000018) |SIEMENS CSA HEADER|
(5528) 0029,0011 Private creator code VR: LO
VM: 1
(5528) (0000021) |SIEMENS MEDCOM HEADER|
(5528) 0029,1008 Modality Image Header Type VR: CS
VM: 1
(5528) (0000009) |NUC FLOOD|
(5528) 0029,1009 Modality Image Header Version VR: LO
VM: 1
(5528) (0000007) |Ver 1.0|
(5528) 0029,1010 Modality Image Header Info VR: OB
VM: 1
(5528) (0032768) <OB, OW, OL, OF data>
(5528) 0029,1131 PMTF Information 1 VR: LO
VM: 1
(5528) (0000013) |25108.0.79319|
(5528) 0029,1132 PMTF Information 2 VR: UL
VM: 1
(5528) (0000004) |262144|
(5528) 0029,1133 PMTF Information 3 VR: UL
VM: 1
(5528) (0000004) |0|
(5528) 0029,1134 PMTF Information 4 VR: CS
VM: 1
(5528) (0000011) |DB TO DICOM|
(5528) 0032,1060 Requested Procedure Description VR: LO
VM: 1
(5528) (0000020) |Myocardial Perfusion|
(5528) 0033,0010 Private creator code VR: LO
VM: 1
(5528) (0000014) |SIEMENS MED NM|
(5528) 0033,1010 COR Data for Detector 1 RETIRED VR: FL
VM: 1-n
(5528) (0000512) |-0.265425|
(5528) |-0.265933|
(5528) |-0.266442|
(5528) |-0.266951|
(5528) |-0.267459|
```


(5528)	-0.267968
(5528)	-0.268477
(5528)	-0.268985
(5528)	-0.269494
(5528)	-0.270003
(5528)	-0.270511
(5528)	-0.27102
(5528)	-0.271529
(5528)	-0.272037
(5528)	-0.272546
(5528)	-0.273055
(5528)	-0.273563
(5528)	-0.274072
(5528)	-0.27458
(5528)	-0.275089
(5528)	-0.275598
(5528)	-0.276106
(5528)	-0.276615
(5528)	-0.277124
(5528)	-0.277632
(5528)	-0.278141
(5528)	-0.27865
(5528)	-0.279158
(5528)	-0.279667
(5528)	-0.280176
(5528)	-0.280684
(5528)	-0.281193
(5528)	-0.281702
(5528)	-0.28221
(5528)	-0.282719
(5528)	-0.283228
(5528)	-0.283736
(5528)	-0.284245
(5528)	-0.284753
(5528)	-0.285262
(5528)	-0.285771
(5528)	-0.286279
(5528)	-0.286788
(5528)	-0.287297
(5528)	-0.287805
(5528)	-0.288314
(5528)	-0.288823
(5528)	-0.289331
(5528)	-0.28984
(5528)	-0.290349
(5528)	-0.290857
(5528)	-0.291366
(5528)	-0.291875
(5528)	-0.292383
(5528)	-0.292892
(5528)	-0.293401
(5528)	-0.293909
(5528)	-0.294418
(5528)	-0.294926
(5528)	-0.295435
(5528)	-0.295944
(5528)	-0.296452

(5528)	-0.296961
(5528)	-0.29747
(5528)	-0.297978
(5528)	-0.298487
(5528)	-0.298996
(5528)	-0.299504
(5528)	-0.300013
(5528)	-0.300522
(5528)	-0.30103
(5528)	-0.301539
(5528)	-0.302048
(5528)	-0.302556
(5528)	-0.303065
(5528)	-0.303574
(5528)	-0.304082
(5528)	-0.304591
(5528)	-0.305099
(5528)	-0.305608
(5528)	-0.306117
(5528)	-0.306625
(5528)	-0.307134
(5528)	-0.307643
(5528)	-0.308151
(5528)	-0.30866
(5528)	-0.309169
(5528)	-0.309677
(5528)	-0.310186
(5528)	-0.310695
(5528)	-0.311203
(5528)	-0.311712
(5528)	-0.312221
(5528)	-0.312729
(5528)	-0.313238
(5528)	-0.313747
(5528)	-0.314255
(5528)	-0.314764
(5528)	-0.315273
(5528)	-0.315781
(5528)	-0.31629
(5528)	-0.316798
(5528)	-0.317307
(5528)	-0.317816
(5528)	-0.318324
(5528)	-0.318833
(5528)	-0.319342
(5528)	-0.31985
(5528)	-0.320359
(5528)	-0.320868
(5528)	-0.321376
(5528)	-0.321885
(5528)	-0.322394
(5528)	-0.322902
(5528)	-0.323411
(5528)	-0.32392
(5528)	-0.324428
(5528)	-0.324937
(5528)	-0.325446

(5528)		-0.325954
(5528)		-0.326463
(5528)		-0.326971
(5528)		-0.32748
(5528)		-0.327989
(5528)		-0.328497
(5528)		-0.329006
(5528)		-0.329515
(5528)		-0.330023
(5528)	0033,1011 COR Data for Detector 2 RETIRED	VR: FL
	VM: 1-n	
(5528)	(0000512)	0.127618
(5528)		0.129177
(5528)		0.130736
(5528)		0.132295
(5528)		0.133854
(5528)		0.135413
(5528)		0.136971
(5528)		0.13853
(5528)		0.140089
(5528)		0.141648
(5528)		0.143207
(5528)		0.144766
(5528)		0.146324
(5528)		0.147883
(5528)		0.149442
(5528)		0.151001
(5528)		0.15256
(5528)		0.154119
(5528)		0.155677
(5528)		0.157236
(5528)		0.158795
(5528)		0.160354
(5528)		0.161913
(5528)		0.163472
(5528)		0.16503
(5528)		0.166589
(5528)		0.168148
(5528)		0.169707
(5528)		0.171266
(5528)		0.172824
(5528)		0.174383
(5528)		0.175942
(5528)		0.177501
(5528)		0.17906
(5528)		0.180619
(5528)		0.182177
(5528)		0.183736
(5528)		0.185295
(5528)		0.186854
(5528)		0.188413
(5528)		0.189972
(5528)		0.19153
(5528)		0.193089
(5528)		0.194648
(5528)		0.196207
(5528)		0.197766

(5528)	0.199325
(5528)	0.200883
(5528)	0.202442
(5528)	0.204001
(5528)	0.20556
(5528)	0.207119
(5528)	0.208678
(5528)	0.210236
(5528)	0.211795
(5528)	0.213354
(5528)	0.214913
(5528)	0.216472
(5528)	0.218031
(5528)	0.219589
(5528)	0.221148
(5528)	0.222707
(5528)	0.224266
(5528)	0.225825
(5528)	0.227384
(5528)	0.228942
(5528)	0.230501
(5528)	0.23206
(5528)	0.233619
(5528)	0.235178
(5528)	0.236737
(5528)	0.238295
(5528)	0.239854
(5528)	0.241413
(5528)	0.242972
(5528)	0.244531
(5528)	0.24609
(5528)	0.247648
(5528)	0.249207
(5528)	0.250766
(5528)	0.252325
(5528)	0.253884
(5528)	0.255442
(5528)	0.257001
(5528)	0.25856
(5528)	0.260119
(5528)	0.261678
(5528)	0.263237
(5528)	0.264795
(5528)	0.266354
(5528)	0.267913
(5528)	0.269472
(5528)	0.271031
(5528)	0.27259
(5528)	0.274148
(5528)	0.275707
(5528)	0.277266
(5528)	0.278825
(5528)	0.280384
(5528)	0.281943
(5528)	0.283501
(5528)	0.28506
(5528)	0.286619

(5528)		0.288178	
(5528)		0.289737	
(5528)		0.291296	
(5528)		0.292854	
(5528)		0.294413	
(5528)		0.295972	
(5528)		0.297531	
(5528)		0.29909	
(5528)		0.300649	
(5528)		0.302207	
(5528)		0.303766	
(5528)		0.305325	
(5528)		0.306884	
(5528)		0.308443	
(5528)		0.310002	
(5528)		0.31156	
(5528)		0.313119	
(5528)		0.314678	
(5528)		0.316237	
(5528)		0.317796	
(5528)		0.319355	
(5528)		0.320913	
(5528)		0.322472	
(5528)		0.324031	
(5528)		0.32559	
(5528)	0033,1014 MHR Y-Shift 1 RETIRED		VR: FL
	VM: 1		
(5528)	(0000004)	-0.053257	
(5528)	0033,1015 MHR Y-Shift 2 RETIRED		VR: FL
	VM: 1		
(5528)	(0000004)	0.053257	
(5528)	0033,1019 NCO Data 2 RETIRED		VR: FL
	VM: 1-n		
(5528)	(0000064)	7.27186	
(5528)		6.95224	
(5528)		6.63619	
(5528)		5.39045	
(5528)		3.70591	
(5528)		2.58014	
(5528)		1.2904	
(5528)		-0.572902	
(5528)		-1.34337	
(5528)		-2.55275	
(5528)		-3.69422	
(5528)		-4.49036	
(5528)		-5.05227	
(5528)		-4.76719	
(5528)		-4.354	
(5528)		-3.7364	
(5528)	0033,101A Non-standard attribute		VR: FL
	VM: 1-n		
(5528)	(0000064)	-4.40101	
(5528)		-6.08608	
(5528)		-7.92421	
(5528)		-9.21146	
(5528)		-9.75772	
(5528)		-10.1813	

```
(5528) | -10.9167 |
(5528) | -11.6636 |
(5528) | -11.7984 |
(5528) | -11.5711 |
(5528) | -10.6393 |
(5528) | -9.48965 |
(5528) | -8.84883 |
(5528) | -7.97996 |
(5528) | -6.76843 |
(5528) | -5.98012 |
(5528) 0033,1020 Bed Correction Angle RETIRED VR: FL
      VM: 1
(5528) (0000004) | 0.0764319 |
(5528) 0033,1021 Gantry Correction Angle RETIRED VR: FL
      VM: 1
(5528) (0000004) | 0.109516 |
(5528) 0033,1022 Bed U/D Correction Data RETIRED VR: SS
      VM: 1-n
(5528) (0000032) | -292 |
(5528) | -239 |
(5528) | -174 |
(5528) | -168 |
(5528) | -210 |
(5528) | -240 |
(5528) | -250 |
(5528) | -259 |
(5528) | -296 |
(5528) | -357 |
(5528) | -462 |
(5528) | -581 |
(5528) | -671 |
(5528) | -800 |
(5528) | -960 |
(5528) | -1108 |
(5528) 0033,1023 Gantry L/R Correction Data RETIRED VR: SS
      VM: 1-n
(5528) (0000032) | -132 |
(5528) | -177 |
(5528) | -229 |
(5528) | -262 |
(5528) | -289 |
(5528) | -349 |
(5528) | -392 |
(5528) | -402 |
(5528) | -485 |
(5528) | -534 |
(5528) | -574 |
(5528) | -615 |
(5528) | -666 |
(5528) | -748 |
(5528) | -807 |
(5528) | -869 |
(5528) 0033,1024 Back Projection Correction Angle head1 VR: FL
      VM: 1
(5528) (0000004) | 0.00306675 |
(5528) 0033,1025 Back Projection Correction Angle head2 VR: FL
      VM: 1
```



```
(5528) |1|
(5528) |1|
(5528) |1|
(5528) |1|
(5528) |1|
(5528) |1|
(5528) 0054,0011 Number of Energy Windows VR: US
      VM: 1
(5528) (0000002) |1|
(5528) 0054,0012 Energy Window Information Sequence VR: SQ
      VM: 1
(5528)
=====
(5528) <Sequence: Item 105 of Message 100 (0054,0012)>
(5528)
=====
(5528) Item ID: 105
(5528) Item Attributes:
(5528) 0054,0013 Energy Window Range Sequence VR:
      SQ VM: 1
(5528)
=====
(5528) <Sequence: Item 106 of Item 105 (0054,0013)>
(5528)
=====
(5528) Item ID: 106
(5528) Item Attributes:
(5528) 0054,0014 Energy Window Lower Limit
      VR: DS VM: 1
(5528) (0000015) |152.33825683594|
(5528) 0054,0015 Energy Window Upper Limit
      VR: DS VM: 1
(5528) (0000015) |177.04174804688|
(5528) ===== End of Item ID: 106 =====
(5528) Item ID: 105 (continued)
(5528)
=====
(5528) <Sequence: Item 107 of Item 105 (0054,0013)>
(5528)
=====
(5528) Item ID: 107
(5528) Item Attributes:
(5528) 0054,0014 Energy Window Lower Limit
      VR: DS VM: 1
(5528) (0000015) |63.665996551514|
(5528) 0054,0015 Energy Window Upper Limit
      VR: DS VM: 1
(5528) (0000015) |77.814002990723|
(5528) ===== End of Item ID: 107 =====
(5528) Item ID: 105 (continued)
(5528) 0054,0018 Energy Window Name VR:
```



```
(5528) 0018,1147 Field of View Shape VR:
      CS VM: 1
(5528) (0000009) |RECTANGLE|
(5528) 0018,1149 Field of View Dimension(s) VR:
      IS VM: 1-2
(5528) (0000006) |421|
(5528) |382|
(5528) 0018,1180 Collimator/grid Name VR:
      SH VM: 1
(5528) (0000009) |@18894496|
(5528) 0018,1181 Collimator Type VR:
      CS VM: 1
(5528) (0000004) |PARA|
(5528) 0018,1182 Focal Distance VR:
      IS VM: 1-2
(5528) (0000001) |0|
(5528) 0018,1183 X Focus Center VR:
      DS VM: 1-2
(5528) (0000002) |0|
(5528) |0|
(5528) 0018,1184 Y Focus Center VR:
      DS VM: 1-2
(5528) (0000002) |0|
(5528) |0|
(5528) 0020,0037 Image Orientation (Patient) VR:
      DS VM: 6
(5528) (0000037) |0.70710670391499|
(5528) |-0.7071068584581|
(5528) |0|
(5528) |0|
(5528) |0|
(5528) |-1|
(5528) 0028,0031 Zoom Factor VR:
      DS VM: 2
(5528) (0000030) |1.4550000429153|
(5528) |1.4550000429153|
(5528) 0054,0200 Start Angle VR:
      DS VM: 1
(5528) (0000014) |225.0030670166|
(5528) ===== End of Item ID: 111 =====
(5528) Message ID: 100 (continued)
(5528)
=====
(5528) <Sequence: Item 112 of Message 100 (0054,0022)>
(5528)
=====
(5528) Item ID: 112
(5528) Item Attributes:
(5528) 0018,1142 Radial Position VR:
      DS VM: 1-n
(5528) (0000239) |216.92630004883|
(5528) |219.03303527832|
(5528) |221.11618041992|
(5528) |229.32727050781|
(5528) |240.43069458008|
```

```
(5528) |247.85102844238|
(5528) |256.35214233398|
(5528) |268.6337890625|
(5528) |273.71221923828|
(5528) |281.68362426758|
(5528) |289.20745849609|
(5528) |294.45504760742|
(5528) |298.15881347656|
(5528) |296.27972412109|
(5528) |293.55627441406|
(5528) |289.48544311523|
(5528) 0018,1145 Center of Rotation Offset VR:
      DS VM: 1
(5528) (0000001) |0|
(5528) 0018,1147 Field of View Shape VR:
      CS VM: 1
(5528) (0000009) |RECTANGLE|
(5528) 0018,1149 Field of View Dimension(s) VR:
      IS VM: 1-2
(5528) (0000006) |421|
(5528) |382|
(5528) 0018,1180 Collimator/grid Name VR:
      SH VM: 1
(5528) (0000009) |@18894497|
(5528) 0018,1181 Collimator Type VR:
      CS VM: 1
(5528) (0000004) |PARA|
(5528) 0018,1182 Focal Distance VR:
      IS VM: 1-2
(5528) (0000001) |0|
(5528) 0018,1183 X Focus Center VR:
      DS VM: 1-2
(5528) (0000002) |0|
(5528) |0|
(5528) 0018,1184 Y Focus Center VR:
      DS VM: 1-2
(5528) (0000002) |0|
(5528) |0|
(5528) 0020,0037 Image Orientation (Patient) VR:
      DS VM: 6
(5528) (0000037) |0.70710682754948|
(5528) |0.70710673482362|
(5528) |0|
(5528) |0|
(5528) |0|
(5528) |-1|
(5528) 0028,0031 Zoom Factor VR:
      DS VM: 2
(5528) (0000030) |1.4550000429153|
(5528) |1.4550000429153|
(5528) 0054,0200 Start Angle VR:
      DS VM: 1
(5528) (0000014) |134.9969329834|
(5528) ===== End of Item ID: 112 =====
(5528) Message ID: 100 (continued)
(5528) 0054,0050 Rotation Vector VR: US
      VM: 1-n
```



```
(5528) |204.14170837402|
(5528) |200.5411529541|
(5528) |197.74940490723|
(5528) |192.90167236328|
(5528) |187.978515625|
(5528) |187.09020996094|
(5528) |188.5888671875|
(5528) |194.73054504395|
(5528) |202.30807495117|
(5528) |206.53193664551|
(5528) |212.25897216797|
(5528) |220.24456787109|
(5528) |225.44058227539|
(5528) 0018,1143 Scan Arc VR:
      DS VM: 1
(5528) (0000002) |90|
(5528) 0018,1144 Angular Step VR:
      DS VM: 1
(5528) (0000005) |5.625|
(5528) 0018,1242 Actual Frame Duration VR:
      IS VM: 1
(5528) (0000005) |60022|
(5528) 0054,0053 Number of Frames in Rotation VR:
      US VM: 1
(5528) (0000002) |16|
(5528) 0054,0200 Start Angle VR:
      DS VM: 1
(5528) (0000014) |225.0030670166|
(5528) ===== End of Item ID: 113 =====
(5528) Message ID: 100 (continued)
(5528) 0054,0090 Angular View Vector VR: US
      VM: 1-n
(5528) (0000064) |1|
(5528) |2|
(5528) |3|
(5528) |4|
(5528) |5|
(5528) |6|
(5528) |7|
(5528) |8|
(5528) |9|
(5528) |10|
(5528) |11|
(5528) |12|
(5528) |13|
(5528) |14|
(5528) |15|
(5528) |16|
(5528) |1|
(5528) |2|
(5528) |3|
(5528) |4|
(5528) |5|
(5528) |6|
(5528) |7|
(5528) |8|
(5528) |9|
```

```
(5528) |10|
(5528) |11|
(5528) |12|
(5528) |13|
(5528) |14|
(5528) |15|
(5528) |16|
(5528) 0054,0202 Type of Detector Motion VR: CS
      VM: 1
(5528) (0000014) |STEP AND SHOOT|
(5528) 0054,0410 Patient Orientation Code Sequence VR: SQ
      VM: 1
(5528)
=====
(5528) <Sequence: Item 101 of Message 100 (0054,0410)>
(5528)
=====
(5528) Item ID: 101
(5528) Item Attributes:
(5528) 0008,0100 Code Value VR:
      SH VM: 1
(5528) (0000007) |F-10450|
(5528) 0008,0102 Coding Scheme Designator VR:
      SH VM: 1
(5528) (0000005) |99SDM|
(5528) 0008,0104 Code Meaning VR:
      LO VM: 1
(5528) (0000009) |Recumbent|
(5528) 0054,0412 Patient Orientation Modifier Code Sequence VR:
      SQ VM: 1
(5528)
=====
(5528) <Sequence: Item 102 of Item 101 (0054,0412)>
(5528)
=====
(5528) Item ID: 102
(5528) Item Attributes:
(5528) 0008,0100 Code Value
      VR: SH VM: 1
(5528) (0000007) |F-10340|
(5528) 0008,0102 Coding Scheme Designator
      VR: SH VM: 1
(5528) (0000005) |99SDM|
(5528) 0008,0104 Code Meaning
      VR: LO VM: 1
(5528) (0000015) |Supine position|
(5528) ===== End of Item ID: 102 =====
(5528) Item ID: 101 (continued)
(5528) ===== End of Item ID: 101 =====
(5528) Message ID: 100 (continued)
(5528) 0054,0414 Patient Gantry Relationship Code Sequence VR: SQ
      VM: 1
(5528)
```

```
=====
(5528) <Sequence: Item 103 of Message 100 (0054,0414)>
(5528)
=====
(5528) Item ID: 103
(5528) Item Attributes:
(5528) 0008,0100 Code Value VR:
      SH VM: 1
(5528) (0000006) |G-5191|
(5528) 0008,0102 Coding Scheme Designator VR:
      SH VM: 1
(5528) (0000005) |99SDM|
(5528) 0008,0104 Code Meaning VR:
      LO VM: 1
(5528) (0000010) |Feet-first|
(5528) ===== End of Item ID: 103 =====
(5528) Message ID: 100 (continued)
(5528) 0055,0010 Private creator code VR: LO
      VM: 1
(5528) (0000014) |SIEMENS MED NM|
(5528) 0055,107E Collimator Thickness RETIRED VR: FL
      VM: 1-n
(5528) (0000008) |25|
(5528) |25|
(5528) 0055,107F Collimator Angular Resolution RETIRED VR: FL
      VM: 1-n
(5528) (0000008) |0.0495146|
(5528) |0.0495146|
(5528) 0055,10C0 Useful Field of View RETIRED VR: SS
      VM: 1-n
(5528) (0000016) |3|
(5528) |0|
(5528) |61|
(5528) |64|
(5528) |3|
(5528) |0|
(5528) |61|
(5528) |64|
(5528) 7FE0,0010 Pixel Data VR: OW
      VM: 1
(5528) <Callback routine>
(5528) 7FE3,0010 Private creator code VR: LO
      VM: 1
(5528) (0000014) |SIEMENS MED NM|
(5528) 7FE3,1014 Minimum Pixel in Frame VR: OW
      VM: 1
(5528) (0000128) <OB, OW, OL, OF data>
(5528) 7FE3,1015 Maximum Pixel in Frame VR: OW
      VM: 1
(5528) (0000128) <OB, OW, OL, OF data>
(5528) 7FE3,1016 Total Counts for each Frame RETIRED VR: OW
      VM: 1
(5528) (0000128) <OB, OW, OL, OF data>
(5528) ===== End of Message ID: 100 =====
```


MESSAGE RECEIVED

```
(5528) 06-26 13:22:54.76 MC3 T7: Message received on context 13
(5528) 06-26 13:22:54.76 MC3 T7: Group 0x0000 Elements:
(5528) 06-26 13:22:54.76 MC3 T7: (0000,0000) Group 0000 Length
(5528) 06-26 13:22:54.76 MC3 T7: (1): 136
(5528) 06-26 13:22:54.76 MC3 T7: (0000,0002) Affected SOP Class UID
(5528) 06-26 13:22:54.76 MC3 T7: (1): 1.2.840.10008.5.1.4.1.1.20
(5528) 06-26 13:22:54.76 MC3 T7: (0000,0100) Command Field
(5528) 06-26 13:22:54.76 MC3 T7: (1): 32769
(5528) 06-26 13:22:54.76 MC3 T7: (0000,0120) Message ID Being Responded To
(5528) 06-26 13:22:54.76 MC3 T7: (1): 0
(5528) 06-26 13:22:54.76 MC3 T7: (0000,0800) Data Set Type
(5528) 06-26 13:22:54.76 MC3 T7: (1): 257
(5528) 06-26 13:22:54.76 MC3 T7: (0000,0900) Status
(5528) 06-26 13:22:54.76 MC3 T7: (1): 0
(5528) 06-26 13:22:54.76 MC3 T7: (0000,1000) Affected SOP Instance UID
(5528) 06-26 13:22:54.76 MC3 T7: (1):
1.3.12.2.1107.5.6.1.6535.32510109062617003339000000002
```

CLEANUP

```
(5556) 06-26 13:22:54.78 MC3 I: CleanupEmptyBlocks: 49 CPlex handled
(5556) 06-26 13:22:54.78 MC3 I: 1 free CPlex found
(5556) 06-26 13:22:54.78 MC3 I: 1792 bytes freed
(5556) 06-26 13:22:54.78 MC3 I: CleanupEmptyBlocks: 47 CPlex handled
(5556) 06-26 13:22:54.78 MC3 I: 39 free CPlex found
(5556) 06-26 13:22:54.78 MC3 I: 160680 bytes freed
(5556) 06-26 13:22:54.78 MC3 I: CleanupEmptyBlocks: 19 CPlex handled
(5556) 06-26 13:22:54.78 MC3 I: 18 free CPlex found
(5556) 06-26 13:22:54.78 MC3 I: 1296 bytes freed
(5556) 06-26 13:22:54.78 MC3 I: CleanupEmptyBlocks: 10 CPlex handled
(5556) 06-26 13:22:54.78 MC3 I: 9 free CPlex found
(5556) 06-26 13:22:54.78 MC3 I: 5472 bytes freed
(5556) 06-26 13:22:54.78 MC3 I: Total bytes freed: 169240
```